JUVÉDERM® VOLUMA® XC

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed physician or properly licensed practitioner.

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

1. DEVICE DESCRIPTION

JUVÉDERM[®] VOLUMA[®] XC is a sterile, biodegradable, non-pyrogenic, viscoelastic, clear, colorless, homogeneous gel implant. It consists of hyaluronic acid (HA) produced by *Streptococcus* species of bacteria, which is crosslinked with BDDE, formulated to a concentration of 20 mg/mL with 0.3% w/w lidocaine in a physiologic buffer.

2. INTENDED USE/INDICATIONS

JUVÉDERM® VOLUMA® XC is indicated for deep (subcutaneous and/or supraperiosteal) injection for cheek augmentation to correct age-related volume deficit in the mid-face and for augmentation of the chin region to improve the chin profile, and for supraperiosteal injection to augment the temple region to improve moderate to severe temple hollowing in adults over the age of 21.

3. CONTRAINDICATIONS

- JUVÉDERM[®] VOLUMA[®] XC is contraindicated for patients with severe allergies manifested by a history of anaphylaxis or history or presence of multiple severe allergies.
- JUVÉDERM[®] VOLUMA[®] XC contains trace amounts of Grampositive bacterial proteins and is contraindicated for patients with a history of allergies to such material.
- JUVÉDERM[®] VOLUMA[®] XC contains lidocaine and is contraindicated for patients with a history of allergies to such material.

4. WARNINGS

- The product must not be injected into blood vessels. Introduction of JUVÉDERM® VOLUMA® XC into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction. Take extra care when injecting soft tissue fillers, for example, after insertion of the needle, and just before injection, the plunger rod can be withdrawn slightly to aspirate and verify the needle is not intravascular, inject the product slowly and apply the least amount of pressure necessary. Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur (see Health Care Professional Instructions #23).
- Product use at specific sites in which an active inflammatory process (skin eruptions such as cysts, pimples, rashes, or hives) or infection is present should be deferred until the underlying process has been controlled.
- Treatment site reactions consist mainly of short-term inflammatory symptoms and generally resolve within 3 days to 4 weeks. Refer to the ADVERSE EVENTS section for details.

5. PRECAUTIONS

- JUVÉDERM[®] VOLUMA[®] XC is packaged for single-patient use. Do not resterilize. Do not use if package is open or damaged.
- In order to minimize the risks of potential complications, this
 product should only be used by health care practitioners
 who have appropriate training, experience, and who are
 knowledgeable about the anatomy at and around the site
 of injection.
- Health care professionals should discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- Based on preclinical studies and a toxicological risk assessment, patients should be limited to 20 mL of any JUVÉDERM[®] injectable gel per 60 kg (132 lbs) body mass per year. The safety of injecting greater amounts has not been established.
- The safety and effectiveness for the treatment of anatomic regions other than the mid-face, chin, pre-jowl sulcus, and temple regions have not been established in controlled clinical studies.
- As with all transcutaneous procedures, dermal filler implantation carries a risk of infection. Standard precautions associated with injectable materials should be followed.
- JUVÉDERM® VOLUMA® XC is to be used as supplied. Modification or use of the product outside the Directions for Use may adversely impact the sterility, homogeneity, and performance of the product.
- The safety for use during pregnancy, in breastfeeding females, and in patients with very thin skin in the mid-face region has not been established.
- The safety has been established for use in patients between 35 and 65 years of age for cheek augmentation, 22 and 80 years of age for chin augmentation, and 32 and 82 years of age for improvement of temple hollowing.
- The safety in patients with known susceptibility to keloid formation, hypertrophic scarring, and pigmentation disorders has not been studied.
- JUVÉDERM[®] VOLUMA[®] XC should be used with caution in patients on immunosuppressive therapy.
- Patients who are using substances that can prolong bleeding (such as aspirin, nonsteroidal anti-inflammatory drugs, and warfarin) may, as with any injection, experience increased bruising or bleeding at treatment sites.
- Patients who experience skin injury near the site of JUVÉDERM[®] VOLUMA[®] XC implantation may be at a higher risk for adverse events.
- Patients may experience late onset nodules with use of dermal fillers, including JUVÉDERM[®] VOLUMA[®] XC. Refer to ADVERSE EVENTS section for details.
- The safety and effectiveness of cannula injection of JUVÉDERM[®] VOLUMA[®] XC has only been clinically evaluated with the TSK STERiGLIDE[®] 25G 1 ½" cannula.
- The safety of JUVÉDERM[®] VOLUMA[®] XC with cannula for cheek augmentation has not been established in patients with Fitzpatrick Skin Types V and VI.
- The safety of JUVÉDERM[®] VOLUMA[®] XC with cannula to improve temple hollowing has not been established.
- After use, treatment syringes and needles may be potential biohazards. Handle and dispose of these items in accordance with accepted medical practice and applicable local, state, and federal requirements.
- JUVÉDERM[®] VOLUMA[®] XC injectable gel is a clear, colorless gel without visible particulates. In the event that the content of a syringe shows signs of separation and/or appears cloudy, do not use the syringe; notify Allergan Product Surveillance at (877) 345-5372.

- JUVÉDERM[®] VOLUMA[®] XC should only be used by health care professionals who have appropriate experience and who are knowledgeable about the anatomy and the product for use in deep (subcutaneous and/or supraperiosteal) injection for cheek and chin augmentation, and supraperiosteal injection to improve moderate to severe temple hollowing.
- Failure to comply with the needle attachment instructions could result in needle disengagement and/or product leakage at the LUER-LOK® and needle hub connection.
- Skin laxity of the chin, neck or jaw could obscure the effects of JUVÉDERM® VOLUMA® XC treatment in the chin region. Therefore, in the chin study, the device was not evaluated in subjects with significant skin laxity of the chin, neck or jaw.
- The effect of JUVÉDERM® VOLUMA® XC injection into the chin on facial hair growth has not been studied.

6. ADVERSE EVENTS

A. Clinical Evaluation of JUVÉDERM[®] VOLUMA[®] XC for **Cheek Augmentation**

In the randomized, controlled clinical trial to evaluate the safety and effectiveness of JUVÉDERM® VOLUMA® XC, there were 238 subjects treated with JUVÉDERM® VOLUMA® XC in the mid-face (zygomaticomalar region, anteromedial cheek, and/ or submalar region, see Figure 1) during the primary phase of the study. Touch-up treatments occurred approximately 30 days after initial injection. After the 6-month blinded "no treatment" control period, control subjects were allowed to receive treatment; 32 control subjects were treated in the study. Preprinted diary forms were used by subjects after treatment to record specific signs and symptoms experienced during each of the first 30 days after initial, touch-up, and repeat treatments in each region of the mid-face. Of the 270 subjects who underwent treatment (from both the treatment and control groups), 265 completed the diary forms. A subset of subjects also underwent repeat treatment following completion of the extended follow-up phase of the study, with 162 subjects completing diary forms after repeat treatment. Subjects were instructed to rate each treatment site response listed on the diary as "Mild (barely noticeable)," "Moderate (uncomfortable)," "Severe (severe discomfort)," or "None."

After initial treatment with JUVÉDERM® VOLUMA® XC, 98% of subjects reported experiencing a local treatment site response. Subjects rated treatment site responses as predominantly mild (21.5%) or moderate (59.2%) in severity with a duration of 2 to 4 weeks. For those treatment site responses evaluated as moderate or severe, the median duration as moderate or severe was 2 days, and the median time to complete resolution was 6 days. Based on data from 167 subjects who received repeat treatment, treatment site responses following repeat treatment were less severe, with reduced incidence and duration compared to initial treatment.

Treatment site responses reported by > 5% of subjects after initial treatments are summarized by severity in Table 1 and by duration in Table 2.

Table 1: Treatment Site Responses by Maximum Severity Occurring in > 5% of Subjects After Initial Treatment for Cheek Augmentation (N = 265)

		Severity ^a		
Treatment Site	Total	Mild	Moderate	Severe
Response	% (n/N ^b)	% (n/N)	% (n/N)	% (n/N)
Any Treatment	98.1%	21.5%	59.2%	19.2%
Site Response	(260/265)	(56/260)	(154/260)	(50/260)
Tenderness	92.1%	46.3%	50.0%	3.7%
	(244/265)	(113/244)	(122/244)	(9/244)
Swelling	85.7%	46.7%	43.6%	9.7%
	(227/265)	(106/227)	(99/227)	(22/227)
Firmness	82.3%	37.6%	54.6%	7.8%
	(218/265)	(82/218)	(119/218)	(17/218)
Lumps/Bumps	81.1%	41.4%	48.8%	9.8%
	(215/265)	(89/215)	(105/215)	(21/215)
Bruising	77.7%	37.4%	51.5%	11.2%
	(206/265)	(77/206)	(106/206)	(23/206)
Pain	66.4%	59.1%	38.6%	2.3%
	(176/265)	(104/176)	(68/176)	(4/176)
Redness	66.0%	60.0%	36.0%	4.0%
	(175/265)	(105/175)	(63/175)	(7/175)
Discoloration	41.1%	62.4%	27.5%	10.1%
	(109/265)	(68/109)	(30/109)	(11/109)
Itching	38.5%	70.6%	18.6%	10.8%
	(102/265)	(72/102)	(19/102)	(11/102)

^a Maximum severity reported in the diary. The denominator for percentages by severity is the number of subjects with the corresponding treatment site response.

^b N denotes number of subjects who recorded responses in the diaries after the initial treatment.

Treatment site responses reported by $\leq 5\%$ of subjects included ache, acne, bulge, bumps, cheek larger upon waking up, dry patch, fine wrinkles, injection/needle marks, numbness, pigmentation from treatment, puffiness, rash, scratch near injection point, soreness, tightness, and yellowness.

Table 2: Duration of Treatment Site Responses After Initial Treatment

			Duration ^a			
Treatment Site	Total	1-3 Days	4-7 Days	8-14 Days	15-30 Days	> 30 Days
Response	% (n/Nʰ)	% (n/N)	% (n/N)	% (n/N)	% (n/N)	% (n/N)
Any Treatment	98.1%	8.1%	22.7%	24.6%	24.6%	20.0%
Site Response	(260/265)	(21/260)	(59/260)	(64/260)	(64/260)	(52/260)
Tenderness	92.1%	29.9%	30.7%	27.9%	8.6%	2.9%
	(244/265)	(73/244)	(75/244)	(68/244)	(21/244)	(7/244)
Swelling	85.7%	41.0%	33.0%	17.6%	5.3%	3.1%
	(227/265)	(93/227)	(75/227)	(40/227)	(12/227)	(7/227)
Firmness	82.3%	26.6%	29.8%	20.2%	11.0%	12.4%
	(218/265)	(58/218)	(65/218)	(44/218)	(24/218)	(27/218)
Lumps/Bumps	81.1%	21.4%	22.3%	22.3%	18.1%	15.8%
	(215/265)	(46/215)	(48/215)	(48/215)	(39/215)	(34/215)
Bruising	77.7%	24.8%	30.6%	29.6%	14.6%	0.5%
	(206/265)	(51/206)	(63/206)	(61/206)	(30/206)	(1/206)
Pain	66.4%	56.3%	31.3%	9.7%	2.8%	0%
	(176/265)	(99/176)	(55/176)	(17/176)	(5/176)	(0/176)
Redness	66.0%	59.4%	28.0%	8.6%	2.3%	1.7%
	(175/265)	(104/175)	(49/175)	(15/175)	(4/175)	(3/175)
Discoloration	41.1%	64.2%	19.3%	6.4%	5.5%	4.6%
	(109/265)	(70/109)	(21/109)	(7/109)	(6/109)	(5/109)
Itching	38.5%	81.4%	16.7%	2.0%	0%	0%
	(102/265)	(83/102)	(17/102)	(2/102)	(0/102)	(0/102)

^a Maximum duration reported in the diary. The denominator for percentages by duration is the number of subjects with the corresponding treatment site response. ^b N denotes number of subjects who recorded responses in the diaries after the initial treatment.

Treatment site responses reported in subject diaries that lasted longer than 30 days were considered adverse events (AEs). AEs were also reported by the Treating Investigator at all follow-up visits where applicable. Table 3 summarizes device/injection-related AEs that occurred with a frequency of > 1%. These adverse events were seen more frequently in subjects that received injection volumes greater than 9 mL and in older subjects (> 60 years). Rarely, adverse events occurred weeks to months after the injection procedure.

Among the 270 treated subjects, 32.6% (88/270) experienced device/injection-related AEs following initial and touch-up treatment, 99% (624/627) of which were reported at a treatment site. The treatment site AEs were evenly divided across the 3 mid-facial regions. Fewer AEs occurred after repeat treatment than after initial/ touch-up treatment.

Table 3: Device/Injection-Related Adverse Events Occurring in > 1% of Treated Subjects for Cheek Augmentation (N = 270)

Adverse Event	Treated Subjects % (n/N)
Treatment site mass	18.9% (51/270)
Treatment site induration	14.1% (38/270)
Treatment site swelling	7.0% (19/270)
Treatment site pain	5.9% (16/270)
Treatment site hematoma	3.7% (10/270)
Treatment site discoloration	2.2% (6/270)
Treatment site erythema	1.9% (5/270)
Treatment site reaction	1.5% (4/270)

Device/injection-related adverse events occurring in $\leq 1\%$ of subjects included injection site hypertrophy (0.7%), nodule (0.7%), inflammation (0.4%), injection site anesthesia (0.4%), injection site dryness (0.4%), injection site erosion (0.4%), mass (0.4%), contusion (0.4%) and syncope (0.4%).

Two subjects (0.7%; 2/270) reported 3 serious adverse events (SAEs) that were considered to be related to the device. Approximately 6 months after treatment, after being scratched near the treated area by a tree branch, one subject experienced inflammation under the left eye. The subject also experienced nodularity in the right cheek approximately 7 months after treatment. The second subject experienced lumps in the cheeks approximately 7 months after treatment. A couple of days before the onset, the subject experienced myofascial pain and body aches. Treatment of the SAEs included topical steroids, oral antibiotics, intralesional steroids, anti-inflammatory medication, and hyaluronidase. All events resolved.

B. 1-Year Post Approval Study of JUVÉDERM[®] VOLUMA[®] XC for Cheek Augmentation

The post-approval study was a statistical evaluation of safety data collected in the JUVÉDERM® VOLUMA® XC pivotal study. Safety data were analyzed from subjects who elected to undergo repeat treatment with JUVÉDERM® VOLUMA® XC as part of the pivotal study. Preprinted diary forms were used by subjects to record specific signs and symptoms experienced during each of the first 30 days after repeat treatment.

Treatment site responses reported by subjects in their diaries after repeat treatment are summarized by severity in Table 4 and by duration in Table 5. The incidence of treatment site responses after repeat treatment was lower than the incidence after initial/touch-up treatment, and treatment site responses were generally less severe and shorter in duration after repeat treatment compared to initial/ touch-up treatment. The majority of treatment site responses after repeat treatment resolved within 2 weeks, while treatment site responses after initial/touch-up treatment typically resolved within 2 to 4 weeks.

Table 4: Treatment Site Responses by Maximum Severity Occurring in > 5% of Subjects After Repeat Treatment for Cheek Augmentation (N = 167)

	Severity ^a								
Treatment Site	Total	Mild	Moderate	Severe					
Response	% (n/N ^b)	% (n/N ^b)	% (n/№)	% (n/N ^b)					
Any Treatment	90.1%	30.8%	54.8%	14.4%					
Site Response	(146/162)	(45/146)	(80/146)	(21/146)					
Tenderness	76.5%	52.4%	42.7%	4.8%					
	(124/162)	(65/124)	(53/124)	(6/124)					
Swelling	67.9%	42.7%	54.5%	2.7%					
	(110/162)	(47/110)	(60/110)	(3/110)					
Firmness	66.0%	40.2%	57.0%	2.8%					
	(107/162)	(43/107)	(61/107)	(3/107)					
Bruising	62.3%	49.5%	37.6%	12.9%					
	(101/162)	(50/101)	(38/101)	(12/101)					
Lumps/Bumps	58.0%	46.8%	47.9%	5.3%					
	(94/162)	(44/94)	(35/94)	(5/94)					
Redness	56.8%	59.8%	38.0%	2.2%					
	(92/162)	(55/92)	(35/92)	(2/92)					
Pain	54.9%	65.2%	30.3%	4.5%					
	(89/162)	(58/89)	(27/89)	(4/89)					
Itching	32.7%	79.2%	20.8%	0%					
	(53/162)	(42/53)	(11/53)	(0/53)					
Discoloration	26.5%	72.1%	27.9%	0%					
	(43/162)	(31/43)	(12/43)	(0/43)					

^a Maximum severity reported in the diary. The denominator for percentages by severity is the number of subjects with corresponding treatment site response.

^b N denotes number of subjects who recorded responses in the diaries after the repeat treatment.

Table 5: Duration of Treatment Site Responses After Repeat Treatment

	Duration ^a								
Treatment Site	Total	1-3 Days	4-7 Days	8-14 Days	15-30 Days	> 30 Days			
Response	% (n/Nʰ)	% (n/Nʰ)	% (n/N ^b)	% (n/Nʰ)	% (n/N⁵)	% (n/N ^b)			
Any Treatment	90.1%	18.5%	30.1%	24.0%	19.2%	8.2%			
Site Response	(146/162)	(27/146)	(44/146)	(35/146)	(28/146)	(12/146)			
Tenderness	76.5%	42.7%	32.3%	13.7%	10.5%	0.8%			
	(124/162)	(53/124)	(40/124)	(17/124)	(13/124)	(1/124)			
Swelling	67.9%	61.8%	23.6%	7.3%	6.4%	0.9%			
	(110/162)	(68/110)	(26/110)	(8/110)	(7/110)	(1/110)			
Firmness	66.0%	26.2%	33.6%	19.6%	14.0%	6.5%			
	(107/162)	(28/107)	(36/107)	(21/107)	(15/107)	(7/107)			
Bruising	62.3%	33.7%	33.7%	23.8%	7.9%	1.0%			
	(101/162)	(34/101)	(34/101)	(24/101)	(8/101)	(1/101)			
Lumps/Bumps	58.0%	37.2%	28.7%	16.0%	10.6%	7.4%			
	(94/162)	(35/94)	(27/94)	(15/94)	(10/94)	(7/94)			
Redness	56.8%	58.7%	29.3%	7.6%	4.3%	0%			
	(92/162)	(54/92)	(27/92)	(7/92)	(4/92)	(0/92)			
Pain	54.9%	65.2%	20.2%	11.2%	3.4%	0%			
	(89/162)	(58/89)	(18/89)	(10/89)	(3/89)	(0/89)			
Itching	32.7%	81.1%	13.2%	5.7%	0%	0%			
	(53/162)	(43/53)	(7/53)	(3/53)	(0/53)	(0/53)			
Discoloration	26.5%	76.7%	7.0%	7.0%	7.0%	2.3%			
	(43/162)	(33/43)	(3/43)	(3/43)	(3/43)	(1/43)			

^a Maximum duration reported in the diary. The denominator for percentages by duration is the number of subjects with corresponding treatment site response.

» N denotes number of subjects who recorded responses in the diaries after the repeat treatment.

Among the 167 subjects who received repeat treatment, 8.4% (14/167) experienced device/injection-related AEs following treatment. All AEs after repeat treatment occurred within 1 month of repeat treatment. The rate of device/injection-related AEs was lower after repeat treatment compared to initial/touch-up treatment. The most common AEs were injection site mass and induration (Table 6).

Table 6: Device/Injection-Related AEs after Repeat Treatment Occurring in > 1% of Treated Subjects for Cheek Augmentation (N = 167)

Adverse Event	Treated Subjects % (n/N)
Injection Site Mass	4.2% (7/167)
Injection Site Induration	4.2% (7/167)
Injection Site Bruising	1.2% (2/167)

All device/injection-related AEs after repeat treatment were mild to moderate, required no action, and resolved without sequelae. Generally, device/injection-related AEs were less severe after repeat treatment compared to initial/touch-up treatment, and most resolved within 3 months. Similar to the initial/touch-up treatment, 3 subjects experienced a device/injection-related AE that lasted more than 180 days, but all resolved without requiring any treatment. Device/injection-related adverse events occurring in \leq 1% of subjects included injection site swelling (0.6%), injection site pain (0.6%), and injection site papule (0.6%).

Of the 121 subjects who completed the 12 months of follow-up after repeat treatment, none experienced any late onset device/ injection-related AEs (those occurring more than 1 month after repeat treatment). There were no device/injection-related serious adverse events after repeat treatment.

C. Other Safety Data JUVÉDERM[®] VOLUMA[®] XC Cannula Study

In the randomized, within-subject controlled clinical trial to evaluate the safety and effectiveness of JUVÉDERM® VOLUMA® XC using cannula in subjects seeking correction of age-related mid-face volume deficit, 60 subjects received treatment using a TSK *STERiGLIDE*® 25G 1 ½" cannula in one cheek and a needle in the other cheek. To achieve optimal correction, the use of a needle was also permitted in the zygomaticomalar region in the cheek randomized to cannula. Preprinted diary forms were used by subjects after treatment to record specific signs and symptoms experienced during each of the first 30 days after treatment. Of the 60 subjects who underwent treatment, 60 completed the diary forms. Subjects were instructed to rate each treatment site response listed on the diary as "Mild (barely noticeable)," "Moderate (uncomfortable)," "Severe (severe discomfort)," or "None."

After treatment with JUVÉDERM® VOLUMA® XC, all subjects reported experiencing a local treatment site response. Subjects rated treatment site responses as being mostly mild or moderate in severity (91.7%), with 60% of subjects having responses resolved within 2 weeks.

Treatment site responses reported by > 5% of subjects after initial treatments are summarized by severity in Table 7 and by duration in Table 8.

AEs were also reported by the Treating Investigator at all follow-up visits, when applicable. Among the 60 mITT subjects, 2 subjects experienced 3 treatment-related AEs (injection site mass on the needle cheek in one subject and injection site plaque on both the needle and cannula cheeks in another subject).

Table 7: Severity of ISRs Occurring in > 5% of Subjects (Safety Population)

Treatment	Treatment Cannula Severity ^a					Needle Severity ^a			
Site	Total	Mild	Moderate	Severe	Total	Mild	Moderate	Severe	
Response	% (n/N⁵)	% (n/N⁵)	% (n/N⁵)	% (n/N⁵)	% (n/N▷)	% (n/N⁵)	% (n/N°)	% (n/N ^b)	
Any Treatment Site Response	100.0% (60/60)	50.0% (30/60)	41.7% (25/60)	8.3% (5/60)	100.0% (60/60)	43.3% (26/60)	45.0% (27/60)	11.7% (7/60)	
Tenderness	91.7%	60.0%	28.3%	3.3%	96.7%	53.3%	38.3%	5.0%	
to touch	(55/60)	(36/60)	(17/60)	(2/60)	(58/60)	(32/60)	(23/60)	(3/60)	
Firmness	83.3%	53.3%	28.3%	1.7%	90.0%	53.3%	31.7%	5.0%	
	(50/60)	(32/60)	(17/60)	(1/60)	(54/60)	(32/60)	(19/60)	(3/60)	
Swelling	81.7%	60.0%	20.0%	1.7%	85.0%	55.0%	28.3%	1.7%	
	(49/60)	(36/60)	(12/60)	(1/60)	(51/60)	(33/60)	(17/60)	(1/60)	
Lumps/	70.0%	51.7%	16.7%	1.7%	83.3%	56.7%	25.0%	1.7%	
Bumps	(42/60)	(31/60)	(10/60)	(1/60)	(50/60)	(34/60)	(15/60)	(1/60)	
Pain after	66.7%	45.0%	18.3%	3.3%	83.3%	56.7%	26.7%	0%	
injection	(40/60)	(27/60)	(11/60)	(2/60)	(50/60)	(34/60)	(16/60)	(0/60)	
Bruising	60.0%	40.0%	16.7%	3.3%	71.7%	41.7%	30.0%	0%	
	(36/60)	(24/60)	(10/60)	(2/60)	(43/60)	(25/60)	(18/60)	(0/60)	
Redness	55.0%	46.7%	8.3%	0%	61.7%	45.0%	16.7%	0%	
	(33/60)	(28/60)	(5/60)	(0/60)	(37/60)	(27/60)	(10/60)	(0/60)	
Discoloration	36.7%	28.3%	8.3%	0%	43.3%	28.3%	15.0%	0%	
	(22/60)	(17/60)	(5/60)	(0/60)	(26/60)	(17/60)	(9/60)	(0/60)	
Itching	18.3%	16.7%	1.7%	0%	20.0%	20.0%	0%	0%	
	(11/60)	(10/60)	(1/60)	(0/60)	(12/60)	(12/60)	(0/60)	(0/60)	

^a Maximum reported severity in the diary.

^b Denominator for percentages is the number of subjects who recorded in the diaries after the treatment.

Table 8: Total Duration of ISRs (Safety Population)

Iu	Table 6. Total Duration of 1945 (Salety Population)									
Cannula Duration ^a					Needle Duration ^a					
Treatment Site Responses	Total % (n/N°)	1-3 Days % (n/№)	4-7 Days % (n/№)	8-14 Days % (n/№)	15-30 Days % (n/N⁵)	Total % (n/N▷)	1-3 Days % (n/N⁵)	4-7 Days % (n/№)	8-14 Days % (n/№)	15-30 Days % (n/№)
Any Treatment Site Response	100.0% (60/60)	16.7% (10/60)	21.7% (13/60)	21.7% (13/60)	40.0% (24/60)	100.0% (60/60)	8.3% (5/60)	25.0% (15/60)	25.0% (15/60)	41.7% (25/60)
Tenderness	91.7%	25.0%	26.7%	23.3%	16.7%	96.7%	23.3%	31.7%	26.7%	15.0%
to touch	(55/60)	(15/60)	(16/60)	(14/60)	(10/60)	(58/60)	(14/60)	(19/60)	(16/60)	(9/60)
Firmness	83.3%	25.0%	21.7%	21.7%	15.0%	90.0%	26.7%	26.7%	21.7%	15.0%
	(50/60)	(15/60)	(13/60)	(13/60)	(9/60)	(54/60)	(16/60)	(16/60)	(13/60)	(9/60)
Swelling	81.7%	33.3%	30.0%	6.7%	11.7%	85.0%	28.3%	35.0%	15.0%	6.7%
	(49/60)	(20/60)	(18/60)	(4/60)	(7/60)	(51/60)	(17/60)	(21/60)	(9/60)	(4/60)
Lumps/	70.0%	20.0%	13.3%	16.7%	20.0%	83.3%	23.3%	23.3%	20.0%	16.7%
Bumps	(42/60)	(12/60)	(8/60)	(10/60)	(12/60)	(50/60)	(14/60)	(14/60)	(12/60)	(10/60)
Pain after injection	66.7%	53.3%	5.0%	5.0%	3.3%	83.3%	53.3%	16.7%	11.7%	1.7%
	(40/60)	(32/60)	(3/60)	(3/60)	(2/60)	(50/60)	(32/60)	(10/60)	(7/60)	(1/60)
Bruising	60.0%	20.0%	18.3%	16.7%	5.0%	71.7%	20.0%	18.3%	20.0%	13.3%
	(36/60)	(12/60)	(11/60)	(10/60)	(3/60)	(43/60)	(12/60)	(11/60)	(12/60)	(8/60)
Redness	55.0%	36.7%	13.3%	5.0%	0%	61.7%	41.7%	16.7%	3.3%	0%
	(33/60)	(22/60)	(8/60)	(3/60)	(0/60)	(37/60)	(25/60)	(10/60)	(2/60)	(0/60)
Discoloration	36.7%	11.7%	11.7%	8.3%	5.0%	43.3%	11.7%	13.3%	11.7%	6.7%
	(22/60)	(7/60)	(7/60)	(5/60)	(3/60)	(26/60)	(7/60)	(8/60)	(7/60)	(4/60)
Itching	18.3%	13.3%	0%	5.0%	0%	20.0%	10.0%	5.0%	3.3%	1.7%
	(11/60)	(8/60)	(0/60)	(3/60)	(0/60)	(12/60)	(6/60)	(3/60)	(2/60)	(1/60)

^a Total calendar days from first to last occurrence of ISR in the diary.

^b Denominator for percentages is the number of subjects who recorded in the diaries after the treatment.

D. JUVÉDERM® VOLUMA® XC Chin Augmentation Study

In a randomized, controlled clinical trial to evaluate the safety and effectiveness of JUVÉDERM® VOLUMA® XC for chin augmentation, there were 144 subjects treated with JUVÉDERM® VOLUMA® XC in the chin area (pogonion, menton, and pre-jowl sulci, see Figure 2) during the primary phase of the study. Touch-up treatments occurred approximately 30 days after initial injection. After the 6-month blinded "no treatment" control period, control subjects could receive treatment; 38 control subjects elected to receive treatment. Subjects were offered repeat treatment 12 months after the last treatment. A total of 74 subjects opted for the repeat treatment. Preprinted diary forms were used by subjects after treatment to record specific signs and symptoms experienced during each of the first 30 days after

initial, touch-up, and repeat treatments. Of the 182 subjects who underwent treatment (from both the treatment and control groups), 181 subjects completed the diary forms, and of the 74 subjects who received repeat treatment, 73 completed the diary forms. Subjects were instructed to rate each treatment site response listed on the diary as "Mild (easily tolerated)," "Moderate (affecting daily activity)," "Severe (unable to do daily activity)," or "None."

After initial treatment with JUVÉDERM® VOLUMA® XC, 92% of subjects reported experiencing a local treatment site response (TSR). Subjects rated TSRs as predominantly mild in severity with a majority (64.7%, 108/167) resolving within 2 weeks. The incidence, severity, and duration of TSRs following repeat treatment were similar to that following initial treatment.

TSRs reported by > 5% of subjects after initial treatments are summarized by severity in Table 9 and by duration in Table 10.

Table 9: Treatment Site Responses by Maximum Severity Occurring in > 5% of Subjects After Initial Treatment for Chin Augmentation (N = 181)

Treatment	Severity ^a						
Site	Total	Mild	Moderate	Severe			
Response	% (n/N⁵)	% (n/N)	% (n/N)	% (n/N)			
Any Treatment Site Response	92.3% (167/181)	44.9% (75/167)	43.1% (72/167)	12.0% (20/167)			
Tenderness	81.8%	56.8%	35.8%	7.4%			
	(148/181)	(84/148)	(53/148)	(11/148)			
Firmness	75.1%	58.8%	36.0%	5.1%			
	(136/181)	(80/136)	(49/136)	(7/136)			
Swelling	68.5%	64.5%	30.6%	4.8%			
	(124/181)	(80/124)	(38/124)	(6/124)			
Pain	63.0%	67.5%	26.3%	6.1%			
	(114/181)	(77/114)	(30/114)	(7/114)			
Lumps/	60.2%	67.9%	25.7%	6.4%			
Bumps	(109/181)	(74/109)	(28/109)	(7/109)			
Bruising	59.1%	59.8%	31.8%	8.4%			
	(107/181)	(64/107)	(34/107)	(9/107)			
Redness	48.6%	69.3%	28.4%	2.3%			
	(88/181)	(61/88)	(25/88)	(2/88)			
Itching	27.6%	86.0%	14.0%	0%			
	(50/181)	(43/50)	(7/50)	(0/50)			
Discoloration	14.9%	74.1%	18.5%	7.4%			
	(27/181)	(20/27)	(5/27)	(2/27)			

^a Maximum severity reported in the diary. The denominator for percentages by severity is the number of subjects with the corresponding treatment site response.

^b N denotes number of subjects who recorded responses in the diaries after the initial treatment.

Table 10: Duration of Treatment Site Responses After InitialTreatment for Chin Augmentation (N = 181)

Treatment	Durationa							
Site	Total	1-3 Days	4-7 Days	8-14 Days	15-30 Days			
Response	% (n/№)	% (n/N)	% (n/N)	% (n/N)	% (n/N)			
Any Treatment Site Response	92.3% (167/181)	13.2% (22/167)	24.0% (40/167)	27.5% (46/167)	35.3% (59/167)			
Tenderness	81.8%	31.8%	43.2%	18.2%	6.8%			
	(148/181)	(47/148)	(64/148)	(27/148)	(10/148)			
Firmness	75.1%	35.3%	27.2%	20.6%	16.9%			
	(136/181)	(48/136)	(37/136)	(28/136)	(23/136)			
Swelling	68.5%	53.2%	31.5%	8.9%	6.5%			
	(124/181)	(66/124)	(39/124)	(11/124)	(8/124)			
Pain	63.0%	69.3%	21.9%	6.1%	2.6%			
	(114/181)	(79/114)	(25/114)	(7/114)	(3/114)			
Lumps/	60.2%	23.9%	21.1%	21.1%	33.9%			
Bumps	(109/181)	(26/109)	(23/109)	(23/109)	(37/109)			
Bruising	59.1%	18.7%	46.7%	31.8%	2.8%			
	(107/181)	(20/107)	(50/107)	(34/107)	(3/107)			
Redness	48.6%	61.4%	22.7%	11.4%	4.5%			
	(88/181)	(54/88)	(20/88)	(10/88)	(4/88)			
Itching	27.6%	70.0%	20.0%	8.0%	2.0%			
	(50/181)	(35/50)	(10/50)	(4/50)	(1/50)			
Discoloration	14.9%	63.0%	22.2%	3.7%	11.1%			
	(27/181)	(17/27)	(6/27)	(1/27)	(3/27)			

^a Maximum duration reported in the diary. The denominator for percentages by duration is the number of subjects with the corresponding treatment site response.

^b N denotes number of subjects who recorded responses in the diaries after the initial treatment.

Overall, 167 treated participants (92.3%) reported at least 1 TSR after initial treatment, 86 (82.7%) reported at least 1 TSR after touch-up treatment, and 55 (75.3%) participants reported at least 1 TSR after repeat treatment.

TSRs reported by $\leq 5\%$ of subjects included pimples, flakiness, numbness, throbbing, tightness, and tingling. TSRs were reported by subjects in their diary for 30 days after treatment. TSRs were considered severe in 12.0% of subjects and lasted 15-30 days in 35.3% of subjects.

Adverse events (AEs) were reported by the TI at all follow-up visits, where applicable. An AE was defined in accordance with ISO 14155 as "any untoward medical occurrence, unintended disease or injury, or untoward clinical signs (including abnormal laboratory findings) in subjects, users, or other persons, whether or not related to the investigational medical device." A treatment-emergent AE (TEAE) was defined as an AE that initially occurred or increased in severity on or after the treatment start date for the treatment group and on or after the randomization date for the control group, and was reported by the treating investigator.

Among the 182 treated subjects, 63 treated participants (34.6%) had 111 TEAEs, and 14 treated participants (7.7%, 14/182) had 20 treatment-related TEAEs following initial and touch-up treatment.

Table 11 summarizes treatment-related TEAEs that occurred with a frequency of > 1%. Treatment-related TEAEs occurring in \leq 1% of subjects included injection site bruising, indentation, induration, inflammation, mass, edema, abscess, cellulitis, gingival pain, and cystic acne, all occurring in 0.5% (1/182) of subjects.

Table 11: Treatment-Related TEAEs Occurring in > 1% of Treated Subjects for Chin Augmentation (N = 182)

Adverse Event	Treated Subjects % (n/N)
Treatment site erythema	1.6% (3/182)
Treatment site pain	1.6% (3/182)

The majority of the treatment-related TEAEs after initial/touch-up treatment were mild or moderate in severity. For initial/touch-up treatment, 2.7% of participants had mild TEAEs, 4.4% had moderate TEAEs, and 1.1% had severe TEAEs.

Fifty percent (7/14) of the participants who experienced treatmentrelated TEAEs resolved within 1 week (Table 13). For initial/touch-up treatment, 3 participants (1.6%) had 4 treatment-related TEAEs that lasted longer than 30 days, including injection site inflammation that lasted 153 days and injection site cellulitis that lasted 36 days, injection site erythema that lasted 264 days, and acne cyst that lasted 134 days.

Fewer AEs occurred after repeat treatment than after initial/touch-up treatment (Table 12). Among the 74 subjects who received repeat treatment, 8 treated participants (10.8%; 8/74) had 12 TEAEs, and 3 treated participants (4.1%; 3/74) had 7 treatment-related TEAEs. The most common TEAE occurring after repeat treatment was injection site mass (2.7%; 2/74). For repeat treatment, 4.1% of participants had mild TEAEs, 1.4% had moderate TEAEs, and 0% had severe TEAEs (Table 13). All TEAEs after repeat treatment did not require any intervention and most resolved within 30 days without sequelae. After repeat treatment, 1 participant (1.4%) had 1 treatment-related TEAE that lasted longer than 30 days: injection site mass that lasted 42 days. There were no serious TEAEs after repeat treatment.

Table 12: Summary of TEAEs after Repeat Treatment (Safety Population)

	Initial Tr	eatment®	Repeat T	reatment	
	Participants (N = 74) Number (%)	Events (N = 26) Number (%)	Participants (N = 74) Number (%)	Events (N = 12) Number (%)	
All TEAEs	21 (28.4%)	26	8 (10.8%)	12	
Treatment-related TEAEs	12 (16.2%)	14	3 (4.1%)	7	
At injection site	10 (13.5%)	11	2 (2.7%)	5	
Not at injection site	2 (2.7%)	3	1 (1.4%)	2	
All SAEs	2 (2.7%)	2	0	0	
Treatment-related SAEs	0	0	0	0	
Discontinued due to TEAE	0	0	0	0	
Deaths	0	0	0	0	

^aAEs with onset within 30 days of initial treatment are included for participants who received repeat treatment to compare the AE rate within the same subset of subjects.

For both initial/touch-up treatments and repeat treatment, most treatment-related TEAEs began within 7 days of treatment. For initial/ touch-up treatment, 1 participant had 3 treatment-related TEAEs that began > 30 days after treatment: injection site edema that began 173 days, 248 days, and 252 days after treatment. These were resolved within 3 days with medication.

There were no treatment related TEAEs that began > 30 days after repeat treatment. All treatment-related TEAEs resolved without sequelae during the study period (12 months of follow-up after initial treatment or 1 month of follow-up after repeat treatment, if applicable). For initial/touch-up treatment, 4 participants had 5 treatment-related TEAEs that required treatment with medication or procedure.

Table 13: Summary of Treatment-related TEAEs for All Treated Participants (Safety Population)

		Number (%)			
	Initial and Touch-up	Treatment	Repeat Ti	reatment	
	Participants	Events	Participants	Events	
	(N = 182)	(N = 20)	(N = 74)	(N = 7)	
Overall Duration	14 (7.7%)	20	3 (4.1%)	7	
≤ 7 days	7 (3.8%)	11	1 (1.4%)	2	
8-14 days	2 (1.1%)	3	0	0	
15-30 days	2 (1.1%)	2	1 (1.4%)	4	
> 30 days	3 (1.6%)	4	1 (1.4%)	1	
Not yet resolved	0	0	0	0	
Time to Onse	t On/After Treatment				
≤ 7 days	12 (6.6%)	15	3 (4.1%)	7	
8-14 days	1 (0.5%)	2	0	0	
15-30 days	0	0	0	0	
> 30 days	1 (0.5%)	3	0	0	
Severity					
Mild	5 (2.7%)	6	3 (4.1%)	6	
Moderate	8 (4.4%)	11	1 (1.4%)	1	
Severe	2 (1.1%)	3	0	0	
Outcome					
Recovered/ Resolved	14 (7.7%)	20	3 (4.1%)	7	
Treatment Re	equired				
No	12 (6.6%)	15	3 (4.1%)	7	
Medication	4 (2.2%)	5	0	0	
Procedure	1 (0.5%)	1	0	0	

Only needle treatment was allowed in the pogonion whereas all participants in the cannula treatment subgroup had some treatment with the needle. Results (Table 14) showed lower incidence of TSRs for injections with cannula than without cannula after each treatment (initial, touch-up, and repeat).

Table 14: Incidence of TSRs After Initial Treatment With and Without Cannula (Safety Population)

TSR	All Treated With Cannula (N = 44) n (%)°	All Treated Without Cannula $(N = 137)$ n $(\%)^{\circ}$
Any TSR	34 (77.3%)	133 (97.1%)
Tenderness to touch	30 (68.2%)	118 (86.1%)
Firmness	28 (63.6%)	108 (78.8%)
Swelling	24 (54.5%)	100 (73.0%)
Bruising	24 (54.5%)	83 (60.6%)
Pain after injection	23 (52.3%)	91 (66.4%)
Lumps/Bumps	21 (47.7%)	88 (64.2%)
Redness	16 (36.4%)	72 (52.6%)
Itching	10 (22.7%)	40 (29.2%)
Discoloration	5 (11.4%)	22 (16.1%)
Other	3 (6.8%)	25 (18.2%)

^aNumber of participants who recorded in the diaries after the treatment.

For the Treated Control group, data after receiving initial treatment at Month 6 are included.

For initial/touch-up treatment with cannula, 2 treated participants (4.3%; 2/46) had 2 treatment-related TEAEs, and without cannula 12 treated participants (8.8%; 12/136) had 18 treatment-related TEAEs (Table 15). For repeat treatment with cannula, there were no treatment related TEAEs; without cannula 3 treated participants (5.4%; 3/56) had 7 treatment-related TEAEs.

Table 15: Comparison of Rate of Treatment-related TEAEs in **Participants Treated With and Without Cannula**

Treatment	AEs in Participants Treated With Cannula and Needle % (n/N)	AEs in Participants Treated Only With Needle % (n/N)
Initial/touch-up treatment	4.3% (2/46)	8.8% (12/136)
Repeat treatment	0 (0/18)	5.4% (3/56)

A total of 11 subjects experienced 14 serious adverse events (SAEs) with onset after the study treatment in the VOLUMA-006 study. One subject (0.5%; 1/182) reported 2 SAEs, injection site inflammation and injection site cellulitis, that were considered to be related to the device (Table 16). These events began 7 days after touch-up treatment and were treated with hyaluronidase, antibiotics, steroid, analgesics/narcotics, sedative, anticoagulant, antacid, electrolyte solutions, antihistamine, and anti-inflammatory medicines. Both events resolved without sequelae, in 36 days for the cellulitis and 153 days for the inflammation, and the participant was discontinued from the study due to these events (the SAE of cellulitis required hospitalization). SAEs that were considered to be not treatment-related were non-cardiac chest pain, invasive papillary breast carcinoma, appendicitis, pneumonia, uterine hemorrhage, keratoacanthoma, squamous cell carcinoma, diverticulitis, atypical pneumonia, intraductal proliferative breast lesion, osteoarthritis, and cholecystitis.

Subjects above the median age (51.5 years and older), experienced more total SAEs (9.9%) than subjects younger than 51.5 years old (2.2%) (Table 27).

	Table 16: Summary of Treatment-Related SAEs					
#	SAE Type	Relationship to Treatment	Duration	Treatme		

#	SAE Type	Relationship to Treatment	Duration	Treatment
1	Injection Site Inflammation	Treatment-related	153 days	Hyaluronidase, Bactrim, Biaxin,
2	Injection Site Cellulitis	Treatment-related	36 days	hydromorphone, diphenhydramine, vancomycin, clindamycin, prednisone, oxycodone, hydrocodone acetaminophen, ibuprofen, heparin, proton pump inhibitors, electrolyte solutions, triamcinolone acetonide sulfamethoxazole, trimethoprim, and valium

Procedural Pain

Participants assessed procedural pain (pain during injection) immediately after completion of each treatment on an 11-point scale ranging from 0 (no pain) to 10 (worst pain imaginable). Pain rated after treatment was minimal, with a mean score of 2.3 for the treatment group at each treatment (initial, touch-up, and repeat) and ranged from 0.0 to 8.0, 0.0 to 7.0, and 0.0 to 6.0 for initial, touch-up, and repeat treatments, respectively. Procedural pain was similar for the treated control participants.

Facial Function Assessments

On the Facial Nerve Grading Scale 2.0, participants were assigned an overall facial function score ranging from I to VI based on the score for individual items (movement in each of 4 facial areas [brow, eye, nasolabial fold, and oral commissure] plus overall facial synkinesis). At baseline and all posttreatment timepoints, over 90% of treated participants were scored as I, indicating the best possible facial function. Only 2 scores of III were given (1 at baseline and 1 at Month 1), both in the treatment group. No scores of IV, V, or VI were given at any point in the study.

Facial Sensation Assessments

Responses for the 2-point discrimination test were the distances for which participants indicated they felt 2 distinct points of pressure at the pogonion and halfway between the pogonion and each pre-jowl sulcus, with possible distances of 1 to 10 mm. Results were similar for baseline and all posttreatment follow-up visits.

The light touch assessment determined the smallest filament number for which participants felt the presence of the filament at the pogonion and halfway between the pogonion and each pre-jowl sulcus, where filaments ranged in diameter from 1.65 to 3.61 mm. The majority of responses at baseline and all posttreatment timepoints occurred with the smallest filament: 1.65 mm.

These facial sensation assessments suggest that treatment did not reduce chin area sensitivity at any timepoint throughout the study.

E. JUVÉDERM[®] VOLUMA[®] XC Temple Hollowing Improvement Study

In a randomized, controlled, multicenter clinical study to evaluate the safety and effectiveness of JUVÉDERM® VOLUMA® XC for correction of temple hollowing, 112 subjects were randomized to treatment group and received injections in the temple area during the primary phase of the study; 58 subjects were randomized to no-treatment control group. Touch-up treatments occurred approximately 30 days after initial injection, if needed, After the 3-month blinded "no treatment" control period, control subjects were offered treatment; 53 control subjects elected to receive treatment. A total of 40 treatment-group subjects opted for the optional maintenance treatment, which was offered to the treatment group 13 months after the last treatment.

Subjects used electronic diaries to record specific signs and symptoms of TSRs experienced during the first 14 days after the initial, touch-up, and maintenance treatments. Subjects in the initial treatment group continued to use the electronic diary on evennumbered days from day 16 to day 30 to record specific signs and symptoms of TSRs. Subjects rated each TSR listed on the diary as "Mild (barely noticeable)," "Moderate (uncomfortable)," "Severe (severe discomfort)," or "None."

After initial treatment with JUVÉDERM® VOLUMA® XC, 59.0% of the subjects (95 of 161) with diary entries reported experiencing at least 1 TSR. Subjects rated TSRs as predominantly mild (70.5%, 67 of 95 reported TSRs) or moderate (26.3%, 25 of 95 reported TSRs) in severity, with a majority of all reported TSRs resolving within 3 days (60%, 57 of 95 reported TSRs). The incidence of TSRs for the touch-up and maintenance treatments was lower than that for initial treatment.

TSRs reported by > 5% of subjects after initial treatment are summarized by severity and duration in Table 17.

Table 17: Treatment Site Responses by Maximum Severity and Duration After Initial Treatment in Temples Occurring in > 5% of **Treated Subjects**

			Severity [®]		Duration			
TSR	Incidence % (n/Nª)	Mild % (n/№)	Moderate % (n/Nª)	Severe % (n/N³)	≤ 3 Days % (n/№)	4-7 Days % (n/№)	8-14 Days % (n/Nª)	15-30 Days % (n/№)
Any TSR	59.0% (95/161)	41.6% (67/161)	15.5% (25/161)	1.9% (3/161)	35.4% (57/161)	12.4% (20/161)	7.5% (12/161)	3.7% (6/161)
Pain After Injection	50.9% (82/161)	37.9% (61/161)	11.2% (18/161)	1.9% (3/161)	37.9% (61/161)	8.7% (14/161)	3.7% (6/161)	0.6% (1/161)
Tenderness to Touch	49.7% (80/161)	42.9% (69/161)	5.6% (9/161)	1.2% (2/161)	34.2% (55/161)	9.9% (16/161)	4.3% (7/161)	1.2% (2/161)
Redness	41.0% (66/161)	34.2% (55/161)	5.6% (9/161)	1.2% (2/161)	34.2% (55/161)	4.3% (7/161)	1.9% (3/161)	0.6% (1/161)
Firmness	37.9% (61/161)	33.5% (54/161)	3.7% (6/161)	0.6% (1/161)	29.2% (47/161)	4.3% (7/161)	2.5% (4/161)	1.9% (3/161)
Swelling	30.4% (49/161)	26.1% (42/161)	4.3% (7/161)	0%	21.1% (34/161)	5.0% (8/161)	3.7% (6/161)	0.6% (1/161)
Lumps/ Bumps	26.1% (42/161)	23.0% (37/161)	3.1% (5/161)	0%	20.5% (33/161)	3.7% (6/161)	1.2% (2/161)	0.6% (1/161)
Bruising	16.1% (26/161)	14.9% (24/161)	1.2% (2/161)	0%	11.2% (18/161)	3.1% (5/161)	1.9% (3/161)	0%
Headaches	9.3% (15/161)	6.8% (11/161)	1.9% (3/161)	0.6% (1/161)	6.2% (10/161)	0%	0%	0%
Discoloration	6.8% (11/161)	5.6% (9/161)	1.2% (2/161)	0%	6.2% (10/161)	0.6% (1/161)	0%	0%

^a N denotes the number of subjects who recorded responses in the diaries after initial treatment.

^b Maximum severity reported in the diary.
 ^c Duration is calculated as total days from first to last date of occurrence.

TSRs reported by \leq 5% of subjects included itching, difficulty chewing, difficulty opening mouth, eye twitching, increased sensation, jaw stiffness, jaw tightness, poor cosmetic result, vasodilation, and wrinkles.

Among the 165 subjects treated with JUVÉDERM® VOLUMA® XC, 17.6% (29 subjects) experienced 51 treatment-related TEAEs following initial and touch-up treatment. All of the treatment-related TEAEs after initial/touch-up treatment were either mild (96.1%, 49 of 51 treatment-related TEAEs) or moderate (3.9%, 2 of 51 treatmentrelated TEAEs) in severity.

Table 18 summarizes treatment-related AEs that occurred with a frequency of > 1%. Treatment-related TEAEs occurring in $\le 1\%$ of subjects included implant site pain, injection site edema, injection site joint discomfort, injection site joint movement impairment, injection site paresthesia, mastication disorder and vasodilation, all occurring in 0.6% (1/165) of subjects.

Table 18: Treatment-Related Adverse Events Occurring in > 1% of Treated Subjects for Improvement of Temple Hollowing (N = 165)

dverse Event	Treated Subjects % (n/N)		
ain in Jaw	6.1% (10/165)		
eadache	4.8% (8/165)		
ijection Site Pain	1.8% (3/165)		
ijection Site Discomfort	1.2% (2/165)		
ijection Site Mass	1.2% (2/165)		
emporomandibular Joint Syndrome	1.2% (2/165)		
rismus	1.2% (2/165)		
rismus	1.2%		

The treatment-related TEAEs occurred within 1 week after treatment for 16.4% (27/165) of subjects and resolved without sequalae within 1 week for 15.2% (25/165) of subjects. There were 7 treatmentrelated TEAEs that began 1 week after treatment affecting 3.0% (5 of 165) of the subjects, which were injection site mass, pain in jaw, headache, and mastication disorder. There were no treatment-related TEAEs that began more than 30 days after treatment.

For initial/touch-up treatment, 4 subjects (2.4%, 4/165) had 6 treatment-related TEAEs that lasted longer than 30 days, including mild trismus that lasted 42 days, mild injection site mass in the left temple that lasted 56 days, mild injection site pain in both temples that lasted 32 days, and mild injection site mass in both temples that lasted 66 days. All the treatment-related TEAEs lasting longer than 30 days resolved without intervention. There were no treatment-related SAEs observed in the study.

Fewer treatment-related TEAEs occurred after maintenance treatment than after initial/touch-up treatment (Table 19). Among the 40 subjects who received maintenance treatment, 12.5% (5/40) experienced treatment-related TEAEs. All TEAEs after maintenance treatment were mild or moderate in severity, did not require any intervention, and most resolved within 30 days without sequelae. There were no serious treatment-related TEAEs after maintenance treatment.

Table 19: Summary of TEAEs for Maintenance Treatment (Safety Population)

(Galety i Opulation)					
	Subjects with TEA Maintenance Tre		Subjects with TEAEs After Initial Treatment		
	Subjects (N = 40) n (%)	Events	Subjects (N = 40) n (%)	Events	
All TEAEs	13 (32.5)	18	14 (35.0)	26	
Treatment-related TEAEs	5 (12.5)	7	6 (15.0)	11	
At Injection Site	1 (2.5)	2	1 (2.5)	1	
Not at Injection Site	4 (10.0)	5	6 (15.0)	10	
All SAEs	4 (10.0)	4	1 (2.5)	2	
Treatment-related SAEs	0 (0)	0	0 (0)	0	
Discontinued due to TEAE	0 (0)	0	0 (0)	0	
Deaths	0 (0)	0	0 (0)	0	
COVID-19 infection TEAEs	3 (7.5)	3	2 (5.0)	2	

Procedural Pain

Subjects assessed procedural pain (pain during injection) immediately after completion of initial, touch-up, and maintenance treatments on an 11-point scale ranging from 0 (no pain) to 10 (worst pain imaginable). Pain was minimal, with mean scores in the treatment group of 1.3 (ranging from 0 to 7) for initial treatment, 1.0 (ranging from 0 to 7) for touch-up treatment, and 1.8 (ranging from 0 to 10) for maintenance treatment. Procedural pain was similar for the control subjects who received delayed treatment.

Facial Function Assessment

Subjects assessed facial function utilizing the Jaw Functional Limitation Scale after the completion of initial, touch-up, and maintenance treatments. Jaw function limitation was minimal, with a median score of zero (no limitation) at baseline and a change from baseline of zero at all postbaseline timepoints in overall score and for the 3 subscale scores for mastication, mobility, and verbal and nonverbal communication.

Vision Assessments

Subjects were assessed on visual acuity and functionality after the completion of initial, touch-up, and maintenance treatments. Snellen visual acuity assessment showed that 11 of 156 (7.1%) subjects had a worsening of their visual acuity while 10 of 156 (6.4%) subjects had an improvement of their visual acuity at Month 1 post initial/touch-up treatments. At Month 1 post maintenance treatment, 5 of 36 (13.9%) subjects had worsening in visual acuity. The visual acuity remained the same for all other subjects post initial/touch-up and maintenance treatments. The TI confirmed that the change in Snellen scores was not clinically significant for any subject eyes that worsened. No eyes showed a \geq 3-line worsening in visual acuity at any assessment.

Confrontational visual fields and ocular motility assessments showed that 100% of eyes were full to confrontation and had full duction and version at all timepoints.

F. Other Safety Data

Postmarket Surveillance

JUVÉDERM® VOLUMA® XC without lidocaine has been marketed outside the US since 2005, and JUVÉDERM® VOLUMA® XC (also known as JUVÉDERM VOLUMA® with lidocaine) has been marketed outside the US since 2009 and in the US since 2013.

The following AEs were received from postmarket surveillance for JUVÉDERM® VOLUMA® XC with and without lidocaine with a frequency of 5 events or more and were not observed during the follow-up in the chin (12 months), temple (13 months), or cheek (24 months) US pivotal clinical studies; this includes reports received globally from all sources including scientific journals and voluntary reports. All AEs obtained through postmarket surveillance are listed in order of number of reports received: edema, inflammatory reaction, non-inflammatory nodule, pain, loss/lack of correction, inflammatory nodule/granuloma, unsatisfactory result, hematoma, allergic reaction, bruising, infection, skin discoloration, device migration, neurological symptoms such as increase or decrease in sensation, vascular occlusion, abscess, anxiety, dermatitis, varied injuries, headache, vision abnormalities, drainage, flu-like symptoms, blister overcorrection, scarring, necrosis, bleeding, malaise, cyst, dry skin, acne, dizziness, herpes, autoimmune disorder exacerbation, angioedema, lymphadenopathy, syncope, nausea, dyspnea, extrusion, telangiectasia, depression, vision loss, calcification, anaphylactic reaction, cardiac complications, biofilm, and deeper wrinkle.

Reported treatments include: antibiotics, steroids, antiseptic creams, hyaluronidase, anti-inflammatories, antihistamines, needle aspiration, eye drops, hyperbaric oxygen treatment, laser treatment, ice, massage, warm compress, analgesics, anti-virals, ultrasound therapy, excision, drainage, surgery, immunosuppressants, petroleum jelly, arnica, anticoagulants, anxiolytics, ACE inhibitors and vasodilators.

Vision abnormalities have been reported following injection of JUVÉDERM® VOLUMA® XC, with and without lidocaine, into the nose, glabella, periorbital area, and/or cheek, with a time to onset

ranging from immediate to 1 week following injection. Reported treatments include anticoagulants, sympathomimetics, steroids, and surgery. Outcomes ranged from resolved to ongoing at the time of last contact. Events requiring medical intervention, and events where resolution information is not available, were reported after injection of JUVÉDERM® VOLUMA® XC with and without lidocaine in the highly vascularized areas of the glabella, nose, and periorbital area, which are outside the device indications for use (see Warnings section).

Delayed-onset inflammation near the site of dermal filler injections is one of the known adverse events associated with dermal fillers. Cases of delayed-onset inflammation have been reported to occur at the dermal filler treatment site following viral or bacterial illnesses or infections, vaccinations, or dental procedures. Typically, the reported inflammation was responsive to treatment or resolved on its own.

Adverse reactions should be reported to Allergan Product Surveillance Department at (877) 345-5372.

7. CLINICAL STUDIES

A. Pivotal Study for JUVÉDERM[®] VOLUMA[®] XC for Cheek Augmentation Pivotal Study Design

A multi-center, single-blind, randomized, no-treatment controlled pivotal clinical study was conducted to evaluate the safety and effectiveness of JUVÉDERM® VOLUMA® XC for cheek augmentation to correct age-related volume deficit in the mid-face. Subjects were randomized to treatment or no-treatment control in a 5.3:1 ratio. Treatment group subjects underwent treatment with JUVÉDERM® VOLUMA® XC at the outset of the study. Up to 2 treatments approximately 1 month apart (initial treatment and up to 1 touch-up treatment) were allowed. The Treating Investigator determined the appropriate volume of JUVÉDERM® VOLUMA® XC to be injected in the 3 sub-regions of the mid-face: zygomaticomalar region, anteromedial cheek region, and submalar region, which are depicted in Figure 1. Treatment of the nasolabial folds and periorbital region was prohibited. The no-treatment control subjects had treatment delayed for 6 months.

Figure 1. Mid-Face Regions Treated

Treated subjects returned for routine safety visits with the Treating Investigator at 1, 3, and 6 months after the last treatment during the primary safety and effectiveness phase. All subjects returned for effectiveness follow-up visits with 2 independent Evaluating Investigators (EI) at 1, 3, and 6 months after the last treatment. Els assessed subjects' overall mid-face volume deficit on the validated 6-point photometric Mid-Face Volume Deficit Scale (MFVDS) as well as volume deficit for each of the 3 facial sub-regions. Els also assessed subjects' improvement on the 5-point Global Aesthetic Improvement Scale (GAIS), the 5-point photometric Nasolabial Fold Photo Severity Scale (NLFSS), and the 11-point Other Aesthetic Features of the Mid-Face guestionnaire. Subjects performed self-assessments on MFVDS, GAIS, NLFSS, treatment goal achievement, satisfaction with mid-facial regions, self-perception of age, look and feel of the face, and satisfaction with facial appearance. Further, 3D facial photography was performed, and volume changes were calculated.

During the extended follow-up period, subjects returned for safety and effectiveness evaluations at quarterly intervals up to 24 months or until any visit at or after Month 12 when the average of the Els' live assessments of the MFVDS returned to, or was worse than, the pre-treatment level. Control subjects followed a similar effectiveness evaluation schedule through Month 6 but were not treated and not required to undergo safety evaluations or self-assessments of effectiveness. After Month 6, control subjects received treatment and followed the same treatment and follow-up schedule as the treatment group. An optional repeat treatment was offered to all subjects after completion of the extended follow-up period, with continued follow-up through 12 months after repeat treatment.

Study Endpoints

The primary effectiveness measure was the average of the 2 blinded Els' live assessments of the subject's overall mid-face volume deficit on the validated 6-point photometric MFVDS. A responder was defined as a subject with \geq 1 grade improvement in the average MFVDS score since baseline.

Effectiveness of JUVÉDERM® VOLUMA® XC was demonstrated if at least 70% of subjects treated with JUVÉDERM® VOLUMA® XC were responders at Month 6, and if the responder rate for the treatment group was statistically superior to that of the no-treatment control group at Month 6.

Secondary measures included the level of improvement on the GAIS and MFVDS assessments for each region of the mid-face as assessed by the blinded Els.

Subject Demographics

A total of 345 subjects were enrolled in the study: 16 were screen failures primarily due to ineligibility, 30 were run-in subjects, and 299 were randomized per protocol, 17 of whom discontinued prior to treatment. Of the remaining 282 subjects, 235 were randomized to the treatment group, and 47 were randomized to the control group. Three-fourths (74.0%, 174/235) of the treatment group completed the extended follow-up period. Sixty-one subjects (26.0%, 61/235) discontinued the study primarily due to loss to follow-up (34.4%, 21/61) or withdrawal of consent (36.1%, 22/61).

At baseline, the majority of subjects in the treatment group (93.6%, 220/235) and all subjects in the control group (100%, 46/46) had moderate, significant, or severe volume deficit (encompassing scores of 2.5 through 5 on the MFVDS scale) in their mid-face according to the average of El assessments.

Subject demographics and pre-treatment characteristics are presented in Table 20.



Characteristic	;	Treatment Group (N = 235) % (n)	Control Group (N = 47) % (n)
Gender	Female	80% (189)	79% (37)
Gender	Male	20% (46)	21% (10)
	Median	56	55
Age (years)	Range (min, max)	(35-65)	(36-65)
	Caucasian	58% (137)	60% (28)
	Hispanic	15% (35)	9% (4)
Race	African-American	19% (44)	26% (12)
	Asian	4% (9)	6% (3)
	Other	4% (10)	0% (0)
	1	3% (6)	4% (2)
	Ш	26% (62)	21% (10)
Fitzpatrick	Ш	29% (67)	23% (11)
Skin Type	IV	18% (43)	30% (14)
	V	19% (44)	19% (9)
	VI	6% (13)	2% (1)

Table 20: Demographics and

Pre-treatment Characteristics (N = 282)

Treatment Characteristics

Multiple injection techniques were used for 95% of subjects, with the most common being tunneling, fanning, and serial puncture. Subjects were injected equally in the 3 facial sub-regions for a total median volume of 2.0 mL for the zygomaticomalar region, 2.0 mL for the anteromedial cheek, and 2.1 mL for the submalar region. The overall total volume used to achieve optimal correction for all 3 sub- regions ranged from 1.2 mL to 13.9 mL, with a median of 6.6 mL. The median volume at initial treatment was 4.8 mL. A touch-up treatment was performed for 82% (195/238) of subjects. The median total volume used for touch-up treatment was 1.9 mL. The median volume injected for repeat treatment was 2.0 mL. The volume of JUVÉDERM[®] VOLUMA[®] XC varied depending on the subject's volume deficit and treatment goal.

Primary Effectiveness Results

JUVÉDERM® VOLUMA® XC provided a clinically and statistically significant improvement in mid-face volume deficit compared to the no-treatment control group. Primary effectiveness was met in that significantly greater than 70% of subjects in the treatment group were responders (85.6% improved by \geq 1 grade compared with their pre-treatment assessment, *p* < 0.0001 against the 70% responder rate threshold), and the responder rate for the treatment group was significantly greater (*p* < 0.0001) than the responder rate for the control group (a difference of 46.7%) at Month 6 (Table 21). JUVÉDERM® VOLUMA® XC was found to be effective in all Fitzpatrick Skin Types, for males and females, and across the studied age range.

Table 21: Effectiveness Summary Responder Rate at 6 Months Based on Evaluating Investigators'

	Responder Rate at Month 6	p-value
Treatment Group	85.6% (178/208)	< 0.0001
Control Group ^a	38.9% (14/36)	
Difference in Responder Rates (Treatment rate - Control rate)	46.7%	< 0.0001

^a Includes 2 subjects who were treated in error.

Secondary Effectiveness Results

The GAIS responder rate for the treatment group was 82.2% (171/208) at Month 6, where the responder rate was the percent of subjects with a score of \geq 1 (improved or much improved) on the GAIS for overall mid-face volume based on EIs' assessments. At Month 6 the MFVDS responder rate for each of the facial sub-regions was above 75%.

Extended Follow-Up

Table 22 shows the mean MFVDS scores during the extended follow-up period (Months 9 to 24). The mean improvement was clinically significant (\geq 1 point), with the majority of subjects demonstrating improvement.

- 86.6% (181/209) at Month 9
- 85.2% (172/203) at Month 12
- 71.5% (128/179) at Month 18
- 67.1% (112/167) at Month 24

Table 22: Mean MFVDS Scores after Initial/Touch-up Treatment

Visit	N	Mean MFVDS Score	Mean Change Since Baseline
Baseline	235	3.3	N/A
Month 9	209	1.7	1.6
Month 12	203	1.8	1.5
Month 18	179	2.1	1.3
Month 24	167	2.2	1.1

Subject Self-Assessments

Subjects performed numerous self-assessments, including satisfaction with facial appearance, self- perception of age, and NLF severity. At each time point, more than three-fourths of the treatment group subjects demonstrated an improvement in the overall satisfaction with facial appearance since baseline. In addition, the majority of treatment group subjects perceived themselves as looking younger than at baseline, from 76.4% at Month 1 to 55.4% at Month 24. Subjects, on average, reported themselves as looking approximately 5 years younger at Month 6 and 3 years younger at Month 24. Lastly, more than half (57%, 236/414) of the treatment group subjects at Month 6 observed a \geq 1-point improvement in their,NLFs.

B. Post-Approval Study for JUVÉDERM[®] VOLUMA[®] XC for Cheek Augmentation

Post-Approval Study Design

The post-approval study was a statistical evaluation of safety data collected from the pivotal study for JUVÉDERM® VOLUMA® XC. The safety data were analyzed for subjects who elected to undergo repeat treatment with JUVÉDERM® VOLUMA®XC as part of the pivotal study. Treated subjects returned for safety visits with the Treating Investigator at 1, 3, 6, 9, and 12 or 12+ months. There were 14 enrolled study sites and 167 subjects from the pivotal study who received repeat treatment and were included in the post-approval study.

The study objective was to compare device/injection-related AEs before and after repeat treatment for subjects who received both initial/touch-up and repeat treatments.

Subject Accountability

After repeat treatment, 93.4% of subjects (156/167) attended the Month 1 visit. Of the 167 subjects, 127 subjects consented to long-term follow-up. At the end of the study, 95.3% of subjects (121/127) completed the follow-up at 12 or 12+ months.

Subject Demographics

Subject demographics (gender, age, race/ethnicity, and Fitzpatrick Skin Type) for the 167 subjects who received repeat treatment are shown in Table 23. The majority of subjects were female and Caucasian, with a median age of 56. Subjects were distributed across all Fitzpatrick Skin Types, with the majority being Fitzpatrick II, III, and IV.

Table 23: Post-Approval Study Subject Demographics

Characteristic		% (n/N) N = 167
Gender	Female	77.8% (130/167)
dender	Male	22.2% (37/167)
	Mean	54.8
	Standard Deviation	6.70
Age (years)	Median	56.0
	Range (Min, Max)	(37, 65)
	Caucasian	64.1% (107/167)
	Hispanic	10.8% (18/167)
Race	African-American	16.2% (27/167)
	Asian	3.6% (6/167)
	Other	5.4% (9/167)
	1	3.6% (6/167)
	Ш	27.5% (46/167)
Fitzpatrick Skin Type	Ш	28.1% (47/167)
	IV	21.0% (35/167)
	V	15.6% (26/167)
	VI	4.2% (7/167)

Study Safety Findings

The primary endpoint was met, demonstrating that the safety profile of JUVÉDERM[®] VOLUMA[®] XC after repeat treatment was not worse than the safety profile after initial/touch-up treatment. The rate of AEs after repeat treatment [8.4% (14/167)] was significantly lower than the rate of AEs after initial/touch-up treatment [33.5% (56/167)].

No new safety concerns were identified after repeat treatment. The types of treatment site responses and AEs observed after repeat treatment were similar to those after initial/touch-up treatment, but were generally less severe after repeat treatment. Most treatment site responses after repeat treatment were mild to moderate and resolved within 2 weeks. Device/injection-related AEs after repeat treatment were all mild to moderate, required no action, and resolved without sequelae, with most resolving within 3 months. After repeat treatment, no subjects experienced any late onset device/injection-related AEs (those occurring more than 1 month after repeat treatment) and no device/injection-related SAEs occurred.

The results of a multivariate analysis demonstrated that device/ injection-related AE rates were different among clinical sites, increased with higher injection volumes, and were higher in females as compared to males.

Study Strengths/Limitations

The strength of the study is that long-term safety data were collected after repeat treatment, with high subject compliance (95.3%) at the end of the study (12/12+ month visit). A limitation of the study was that the length of follow-up after repeat treatment was approximately one year shorter than after initial/touch-up treatment.

C. JUVÉDERM[®] VOLUMA[®] XC Cannula Study

Study Design

A multi-center, evaluator-blinded, randomized, within-subject controlled clinical study was conducted to evaluate the safety and effectiveness of JUVÉDERM® VOLUMA® XC with cannula to correct age-related mid-face volume deficit. Subjects were randomized to undergo treatment with the TSK *STERiGLIDE*® 25G 1 ½" cannula in one cheek and a needle in the other cheek. The use of a needle was also permitted in the zygomaticomalar region in the cheek randomized to cannula to achieve optimal correction. At the outset of the study, 60 enrolled subjects underwent treatment with JUVÉDERM® VOLUMA® XC.

Treated subjects returned for routine safety visits with the TI at 1 and 3 months after the treatment. At these visits, blinded Els assessed subjects' overall mid-face volume deficit on the validated 6-point photometric MFVDS. Subjects performed self-assessments on the Satisfaction with Cheeks module of the FACE-Q Questionnaire.

Study Endpoints

The primary effectiveness measure was the blinded El's assessment at Month 1 of the subject's volume deficit for each cheek on the validated 6-point photometric MFVDS. The primary effectiveness endpoint was to demonstrate non-inferiority of JUVÉDERM® VOLUMA® XC administered via a TSK *STERiGLIDE*® 25G 1 ½" cannula versus needle.

Secondary measures included EI-assessed overall MFVDS responder rates and subject-assessed mean overall satisfaction scores on the validated *Satisfaction with Cheeks* module of the FACE-Q Questionnaire at Month 1. A responder was defined as a subject with \geq 1 grade improvement in the MFVDS score since baseline.

Subject Demographics

A total of 60 subjects received treatment. At baseline, all subjects had moderate, significant, or severe volume deficit (encompassing scores of 3 through 5 on the MFVDS scale) in their mid-face according to the El's assessments. Subject demographics and pre-treatment characteristics are presented in Table 24.

Table 24: Subject Demographics		
Parameter Total (N =		
Gender		
Female	49 (81.7%)	
Male	11 (18.3%)	
Age		
Median (SD)	54.9 (6.41)	
Range (Min, Max)	37, 65	
Fitzpatrick Skin Type		
1	2 (3.3%)	
	15 (25.0%)	
Ш	31 (51.7%)	
IV	11 (18.3%)	
٧	0	
VI	1 (1.7%)	

Table 24: Subject Demographics

Primary Effectiveness Results

This study demonstrated non-inferiority of JUVÉDERM[®] VOLUMA[®] XC administered via a TSK *STERiGLIDE*[®] 25G 1 ½" cannula versus needle. Based on the blinded El assessments, the 95% confidence interval for the difference in mean change in MFVDS from baseline to Month 1 (cannula minus needle) was 0.1. The upper confidence limit was < 0.5 and therefore, statistical non-inferiority was concluded.

Secondary Effectiveness Results

The EI-assessed MFVDS responder rates at Month 1 were 93.3% for cannula-treated and 95.0% for needle-treated cheeks, with a paired difference of -1.7 (95% CI: -7.31 to 3.98). The mean score at baseline on the *Satisfaction with Cheeks* module of the FACE-Q questionnaire was 32.1 and at Month 1 increased by 55.5 points; over 85% of subjects had improved satisfaction with the attractiveness, youthful fullness, and contour of their cheeks.

Other Effectiveness Results

The El assessed MFVDS responder rates at Month 3 were 93.3% for cannula-treated and 98.3% for needle-treated cheeks, with a paired difference of -5.0 (95% CI: -10.51 to 0.51).

D. Pivotal Study for JUVÉDERM[®] VOLUMA[®] XC for Chin Augmentation

Study Design

A multi-center, single-blind, randomized, no-treatment controlled pivotal clinical study was conducted to evaluate the safety and effectiveness of JUVÉDERM® VOLUMA® XC for chin augmentation. Subjects were randomized to treatment or no-treatment control in a 3:1 ratio. Treatment group subjects underwent treatment with JUVÉDERM® VOLUMA® XC at the outset of the study. The Treating Investigator (TI) determined the appropriate volume of JUVÉDERM[®] VOLUMA[®] XC to be injected in the chin area (did not exceed 4 mL for initial and touch-up treatment combined and another 4 mL for repeat treatment): pogonion, mentum, and pre-jowl sulci, as depicted in Figure 2. Injection in the pogonion was only permitted with a 27G ¹/₂" needle; a TSK 25 G 1 1/2" *STERiGLIDE*[®] cannula was permitted for injection in the other treatment areas. The no-treatment control subjects had treatment delayed for 6 months.

Figure 2. Chin Area Treated



Up to 2 treatments approximately 1 month apart (initial treatment and up to 1 touch-up treatment) were allowed. All subjects returned for routine safety and effectiveness follow-up visits at 1, 3, and 6 months after the last treatment during the primary safety and effectiveness phase. During the extended follow-up period, treatment group subjects returned for safety and effectiveness evaluations at 9 and 12 months after last treatment. An optional repeat treatment was offered to all treatment group subjects after completion of the extended follow-up period, with 1 month of follow up after repeat treatment. Control subjects followed a similar effectiveness evaluation schedule through Month 6. After Month 6, control subjects received treatment and were followed for an additional 6 months with the same treatment and follow-up schedule as the treatment group.

Pre- and post-procedure, the objective parameters measured during the study included the evaluating investigators' (Els') assessment of subjects' overall chin volume deficit live and via 2D profile images of the left side of the chin, which were rendered by image analysis software from 3D photos, using the validated 5-point photonumeric Allergan Chin Retrusion Scale (ACRS). Els also assessed subjects' improvement on the 5-point Global Aesthetic Improvement Scale (GAIS). Subjects performed self-assessments on the GAIS, the *Satisfaction with Chin* module of the validated FACE-Q questionnaire, the Psychological Well-Being module of the validated FACE-Q questionnaire, and the natural look and natural feel of the chin area on an 11-point scale. Further, 3D facial photography was performed to quantify volume changes. In addition to subject diaries and TI assessment of AEs, safety of the treatment was assessed via facial sensation and facial function assessments conducted by the El.

Sensation in the chin area was assessed using a 2-point discrimination test and a light touch test. Two-point discrimination was tested at 3 locations on the chin (the pogonion and halfway between the pogonion and each prejowl sulcus). Using the Dellon Disk-Criminator, the El lightly touched the set of prongs to each location on the chin, while the subject reports whether s/he felt "1" or "2" objects touching his/her skin.

The light touch test was performed at the same 3 locations on the chin. The El pressed Semmes-Weinstein monofilaments of different diameters against the subject's skin and record the smallest filament size that elicits a response at each assessed location.

The EI tested facial function using the Facial Nerve Grading Scale 2.0 (FNGS 2.0). The EI assessed the subject's face at rest and then requested that the subject make a standardized series of facial movements while the EI rated the movement in each of 4 facial areas: brow, eye, nasolabial fold, and oral commissure. A score was assigned to each facial area, and a score for synkinesis is attributed across the entire face.

Study Endpoints

With regards to safety, preprinted diary forms were used by subjects after treatment to record specific signs and symptoms experienced during each of the first 30 days after initial, touch-up, and repeat treatments. Subjects were instructed to rate each treatment site response listed on the diary as "Mild (easily tolerated)," "Moderate (affecting daily activity)," "Severe (unable to do daily activity)," or "None." Adverse Events were reported by the TI at all follow-up visits where applicable. With regards to effectiveness, the primary effectiveness measure was the single blinded EI's assessment of the subject's chin volume deficit in 2D images using the validated 5-point photonumeric ACRS (Table 25, Figure 3).

The ACRS scale was validated in a 61-subject study where three reviewers were shown photographs of the subjects at two different time points. The average weighted kappa for the intra-rater agreement was 0.87, meaning the reviewers' evaluations of the same subjects were consistent between the two time points. The weighted kappa for the inter-rater agreement for two of the three reviewers was 0.84. The agreement between those two reviewers and the third reviewer was lower than 0.6 (0.59 and 0.53), but the balance of evidence suggested that the scale could be used consistently.

Secondary measures included the statistical superiority, at Month 6 compared to baseline, of the mean overall score on the *Satisfaction with Chin* module of the validated FACE-Q questionnaire (0 to 100, where higher scores reflect a better outcome) as assessed by the subjects, and the level of improvement on the GAIS as assessed by the blinded EIs and the subjects. Other effectiveness endpoints included the responder rate and ACRS score based on the EI's live assessment at baseline and Month 6.

With regards to success/failure criteria, a responder was defined as a subject with \ge 1-point improvement in the ACRS score compared to the baseline score. Effectiveness of JUVÉDERM® VOLUMA® XC was demonstrated if at least 50% of subjects treated with JUVÉDERM® VOLUMA® XC were responders (\ge 1-point improvement) at Month 6, and if the responder rate for the treatment group was statistically superior to that of the no-treatment control group at Month 6.

	Table 20. Allergan onin neurasion ooale		
Score	Grade	Description	
0	None	No chin retrusion; Chin midpoint* at or in front of the lower vermilion border vertical line	
1	Minimal	Minimal chin retrusion; Chin midpoint* is between the labiomental sulcus vertical line and lower vermilion border vertical line	
2	Moderate	Moderate chin retrusion; Chin midpoint* at labiomental sulcus vertical line	
3	Severe	Severe chin retrusion; Chin midpoint* slightly behind labiomental sulcus vertical line	
4	Extreme	Extreme chin retrusion; Chin midpoint* significantly behind labiomental sulcus vertical line	

*Chin midpoint: the midpoint between the labiomental sulcus and the inferior point of the chin.

Figure 3: Lines Used in the ACRS Descriptors



Subject Demographics

Data from all 221 enrolled subjects were available for analysis (Table 26). Of the 221 subjects, 29 were screen failures primarily due to ineligibility, and 192 were randomized per protocol, with 144 in the treatment group and 48 in the control group. Of the 192 randomized subjects, 169 (88.0%; 128 treatment and 41 control) completed the Month 6 primary endpoint visit, and 38 of the 48 control group subjects (79.2%) opted to receive study treatment after the completion of the 6-month control period. A total of 167 (87.0%; 127 treatment and 40 control) completed the study.

At baseline, 7.8% (15/192) of subjects had mild, 40.6% (78/192) had moderate, 43.2% (83/192) had severe, and 8.3% (16/192) had extreme chin volume deficit based on El photo assessments on the ACRS.

Table 2	6: Partici	ipant Dis	position
---------	------------	-----------	----------

Population	Nu		
	Treatment Group	Control Group	Total
Enrolled	N/#	Ň.	221
Screen failures	N/A	N N	29
Randomized	144	48	192
Modified Intent-to-Treat (mITT) Population	144	48	192
Number of Participants Treated	144	38	182
Completed Month 6 visit (primary endpoint)	128	41	169
Completed Month 12 visit	128	N/A	128
Received Repeat Treatment at Month 12	74	N/A	74
Completed Study	127	40	167
Discontinued from the Study	17	8	25
Adverse Event*	1	2	3
Lost to follow-up	10	1	11
Personal Reasons	6	4	10
Protocol Deviation**	0	1	1
Per-Protocol (PP) Population	139	48	187
Safety Population	144	48	192

* The one treatment group participant was discontinued from the study due to a treatment-related serious adverse event, as detailed in Section 6D. The two control group participants discontinued due to adverse events prior to receiving any study treatment.

** The one control group participant discontinued from the study due to protocol deviation of unable to comply with the study visit schedule.

The demographics of the study population are typical for a study performed in the US. Subject demographics and pre-treatment characteristics are presented in Table 27.

 Table 27: Demographics and Pre-treatment Characteristics (N = 192)

Characteristic	;	Treatment Group (N = 144) % (n)	Control Group (N = 48) % (n)
Gender	Female	90% (129)	85% (41)
Gender	Male	10% (15)	15% (7)
	Median	51.5	52.5
Age (years)	Range (min, max)	(23-80)	(22-72)
	White	58% (137)	60% (28)
	Black or African-American	15% (35)	9% (4)
Race	Asian	19% (44)	26% (12)
	American-Indian or Alaska Native	4% (9)	6% (3)
	Multiple	4% (10)	0% (0)
	Hispanic or Latino	19% (27)	13% (6)
Ethnicity	Not Hispanic or Latino	81% (117)	87% (42)
	1	5% (7)	2% (1)
	Ш	31% (44)	31% (15)
Fitzpatrick Skin Type	Ш	37% (54)	33% (16)
	IV	15% (22)	17% (8)
	٧	7% (10)	13% (6)
	VI	5% (7)	4% (2)

Treatment Characteristics

The most common injection techniques at any treatment were bolus and serial puncture. At initial treatment, 99.3% of treatment group subjects were treated in the pogonion, 77.8% in the menton, and 87.5% in the pre-jowl sulci. Needles were used for 100% of subjects, and cannulas were used for 25.0% of subjects at initial treatment. The median total volume used to achieve optimal correction was 2.4 mL (range, 0.7-4.0 mL), with 1.0 mL in the pogonion, 0.5 mL in the menton, and 1.0 mL in the pre-jowl sulci (right and left combined). The median volume at initial treatment was 2.0 mL. A touch-up treatment was performed for 45.8% (88/192) of subjects with a median total volume of 1.0 mL. The repeat treatment was performed for 45.8% (88/192) of subjects and the median volume injected for repeat treatment was 2.0 mL. The volume of JUVÉDERM® VOLUMA® XC varied depending on the subject's chin volume deficit and treatment goal.

Primary Effectiveness Results

JUVÉDERM[®] VOLUMA[®] XC provided a clinically and statistically significant improvement in chin volume deficit compared to the no-treatment control group. The analysis of primary effectiveness was based on the 126 treatment group and 40 control group evaluable subjects at the 6-month time point. The primary effectiveness endpoint was met in that greater than 50% of subjects in the treatment group were responders (56.3% improved by \geq 1 point compared with their pre-treatment assessment), and the responder rate for the treatment group was significantly greater (p = 0.0019) than the responder rate for the control group (a difference of 28.8%) at Month 6 (Table 28). The primary effectiveness endpoint was below 50% in the following subgroups: older subjects (aged 51.5 years and older), darker skin types (FST V/VI), and males. However, the satisfaction rates were high and the secondary endpoints were met for these subgroups.

Table 28: Effectiveness Summary Responder Rate at 6 Months Based on Evaluating Investigators' Assessments of Images

	Responder Rate at Month 6	p-value
Treatment Group	56.3% (71/126)	N/A
Control Group	27.5% (11/40)	N/A
Difference in Responder Rates (Treatment rate - Control rate)	28.8%	0.0019

The responder rate at the 12-month follow-up visit was 57.6% (72/125) and 73.9% (51/69) at 1 month after repeat treatment based on the blinded El's assessment of 2D images.

The responder rate based on the live assessment of subjects at Month 6 was 91.8% (89/97) for the treatment group and 23.3% (7/30) for the control group.

Secondary effectiveness results: The FACE-Q Satisfaction with *Chin* overall mean score was 34.9 at baseline and improved to 71.3 at Month 6 with the improvement being statistically significant (p < 0.0001) (Figure 4). Most of the subjects (91.8%) reported satisfaction with their chin 1 month after treatment. Among other questions, this FACE-Q questionnaire included questions on satisfaction with chin look in profile view and width of the chin. At Month 1, 88.8% of treatment group subjects were satisfied with how their chin looks in profile view and 95.5% were satisfied with the width of their chin.



Figure 4: FACE-Q Satisfaction with Chin Mean Scores by Visit

Baseline Month 1 Month 3 Month 6 Month 9 Month 12 M1 After Repeat

Visit

Treatment (Tx) Control (Ctrl)

The EI and subject GAIS responder rates at Month 6 for the treatment group were 91.2% (114/125) and 87.3% (110/126), respectively, where the responder rate was the percent of subjects with a score of improved or much improved on the GAIS. The EI GAIS responder rate at Month 6 for the untreated control group was 19.5% (8/41) for EI.

An independent, blinded assessment was conducted on full-face 3-dimensional (3D) images collected at randomization (baseline) and at follow-up visits, including the primary timepoint (Month 6). Three independent raters used the ACRS to assess the severity of chin retrusion in each 3D image. At Month 6 the mean change in ACRS score for the treatment group was statistically superior to that for the untreated control group (p < 0.0001). However, at Month 6 the ACRS responder rate for the treatment group was less than 50%, though it was greater than that for the untreated control group (43.0% versus 12.5%).

Subgroup Analyses

The following characteristics were evaluated for potential association with outcomes: Fitzpatrick Skin Type (FST) (Table 29), age (Table 30), gender (Table 31), baseline ACRS, injection volume, cannula usage (Table 14 and 15), and investigational site.

Table 29: Effectiveness and Safety Results at 6 Months by FST Subgroups

		Fitzpa	atrick Skin Type Subgroup		
Assessment Group		I/II	III/IV	v/vi	
EFFECTIVENESS*					
2D ACRS	Treatment	63.0% (29/46)	57.6% (38/66)	28.6% (4/14)	
Responder Rate, % (n/N)	Control	12.5% (2/16)	33.3% (6/18)	50.0% (3/6)	
EI GAIS	Treatment	93.5% (43/46)	89.4% (59/66)	92.3% (12/13)	
Responder Rate, % (n/N)	Control	18.8% (3/16)	15.8% (3/19)	33.3% (2/6)	
Subject GAIS Responder Rate, % (n/N)	Treatment	87.0% (40/46)	89.4% (59/66)	78.6% (11/14)	
FACE-Q Satisfaction	Treatment	71.0 (46)	71.7 (66)	70.0 (14)	
with Chin Mean Score (n)	Control	30.9 (16)	30.5 (19)	41.0 (6)	
FACE-Q Satisfaction	Treatment	35.0 (46)	35.8 (66)	36.2 (14)	
with Chin Mean Change from Baseline (n)	Control	-2.7 (16)	-2.9 (18)	-6.2 (6)	
Live ACRS	Treatment	97.1% (33/34)	90.7% (49/54)	77.8% (7/9)	
Responder Rate, - % (n/N)	Control	11.1% (1/9)	27.8% (5/18)	33.3% (1/3)	
Mean Change in Volume in cc	Treatment	2.0 (46)	2.9 (66)	1.6 (14)	
using 3D Image Analysis (n)**	Control	-0.03 (16)	0.53 (18)	-0.6 (6)	
		SAFETY*			
Total TEAEs		33.0% (22/66)	36.6% (34/93)	30.4% (7/23)	
Treatment- related TEAEs		7.6% (5/66)	6.5% (6/93)	13.0% (3/23)	
All SAEs	Treatment	7.6% (5/66)	4.3% (4/93)	8.7% (2/23)	
Treatment Site Responses after Initial Treatment		98.5% (65/66)	90.2% (83/92)	82.6% (19/23)	

* The N for the effectiveness data is only the treatment group and the N for the safety data includes all treated subjects.

** Median injection volume for the treatment group was 2.0 mL, 2.5 mL, and 2.0 mL in FST I/II, III/IV, and V/VI groups, respectively.

For subjects with darker skin (FST V/VI) (Table 29), the device did not meet the primary effectiveness endpoint (28.6% responder rate in the treatment group) and performed worse than the no-treatment control (50.0% responder rate, where a responder is a subject with at least a 1-point improvement in the 2D ACRS score at 6 months from baseline).

Assessment	Group	Age Subgroup		
Assessment	Group	< 51.5 years	≥ 51.5 years	
	EFFECTIVENESS	*	• •	
2D ACDC Deenender Date 9/ (n/N)	Treatment	67.8% (40/59)	46.3% (31/67)	
2D ACRS Responder Rate, % (n/N)	Control	26.3% (5/19)	28.6% (6/21)	
FLOAIC Descender Date 0/ (a/h)	Treatment	91.5% (54/59)	90.9% (60/66)	
El GAIS Responder Rate, % (n/N)	Control	42.1% (8/19)	0 (0/22)	
Subject GAIS Responder Rate, % (n/N)	Treatment	93.2% (55/59)	82.1% (55/67)	
FACE-Q Satisfaction with Chin	Treatment	72.1 (59)	70.5 (67)	
Mean Score (n)	Control	32.7 (19)	31.7 (22)	
FACE-Q Satisfaction with Chin	Treatment	37.4 (59)	33.9 (67)	
Mean Change from Baseline (n)	Control	-4.2 (18)	-2.6 (22)	
Live ACDC Deepender Date (// /n/N)	Treatment	97.9% (46/47)	86.0% (43/50)	
Live ACRS Responder Rate, % (n/N)	Control	30.8% (4/13)	17.6% (3/17)	
Mean Change in Volume in cc using	Treatment	2.8 (59)	2.1 (67)	
3D Image Analysis (n)**	Control	0.02 (19)	0.25 (21)	
	SAFETY*		• •	
Total TEAEs	Treatment	29.7% (27/91)	39.6% (36/91)	
Treatment-related TEAEs	Treatment	7.7% (7/91)	7.7% (7/91)	
All SAEs	Treatment	2.2% (2/91)	9.9% (9/91)	
Treatment Site Responses after Initial Treatment	Treatment	97.8% (88/90)	86.8% (79/91)	

Table 30: Effectiveness and Safety Results at 6 Months by Age Subgroups

* The N for the effectiveness data is only the treatment group and the N for the safety data includes all treated subjects.

** Median injection volume for the treatment group was 2.5 mL and 2.1 mL in < 51.5 and \geq 51.5, respectively.

The median age of subjects was 51.5 years. For subjects aged 51.5 years and older (Table 30), the device did not meet the primary effectiveness endpoint (46.3% responder rate in the treatment group). The responder rate for subjects aged 51.5 years or older is statistically higher than the no-treatment control group (28.6%).

A	0	Gender Subgroup		
Assessment	Group	Female	Male	
	EFFECTIVENESS	*		
2D ACRS Responder Rate, % (n/N)	Treatment	57.7% (64/111)	46.7% (7/15)	
2D ACRS Responder Rate, % (II/N)	Control	32.4% (11/34)	0% (0/6)	
El GAIS Responder Rate, % (n/N)	Treatment	91.9 (102/111)	85.7% (12/14)	
El GAIS Responder rate, % (II/N)	Control	20.0% (7/35)	16.7% (1/6)	
Subject GAIS Responder Rate, % (n/N)	Treatment	86.5% (96/111)	93.3% (14/15)	
FACE-Q Satisfaction with Chin	Treatment	71.3 (111)	71.0 (15)	
Mean Score (n)	Control	31.5 (34)	35.0 (6)	
FACE-Q Satisfaction with Chin	Treatment	35.0 (111)	39.5 (15)	
Mean Change from Baseline (n)	Control	-3.9 (34)	0.2 (6)	
Live ACRS Responder Rate, % (n/N)	Treatment	92.9% (79/85)	83.3% (10/12)	
Live Acho nesponder hate, % (II/N)	Control	26.9% (7/26)	0 (0/4)	
Mean Change in Volume in cc using	Treatment	2.44 (111)	2.39 (15)	
3D Image Analysis (n)**	Control	0.01 (34)	0.85 (6)	
	SAFETY*			
Total TEAEs	Treatment	35.8% (58/162)	25.0% (5/20)	
Treatment-related TEAEs	Treatment	8.0% (13/162)	5.0% (1/20)	
All SAEs	Treatment	6.8% (11/162)	0 (0/20)	
Treatment Site Responses after Initial Treatment	Treatment	93.8% (151/161)	80.0% (16/20)	

Table 31: Effectiveness and Safety Results at 6 Months by Gender Subgroups

* The N for the effectiveness data is only the treatment group and the N for the safety data includes all treated subjects.

** Median injection volume for the treatment group was 2.0 mL and 3.7 mL in females and males, respectively. For male subjects (Table 31), the device did not meet the primary effectiveness endpoint (46.7% responder rate in the treatment group). There were no responders in the male control subjects.

The effectiveness of JUVÉDERM[®] VOLUMA[®] XC for chin augmentation in the VOLUMA-006 study was similar for subjects with:

- Moderate or severe chin retrusion
- Treatment with or without cannula

By investigational site, the responder rate for the treatment group based on photo assessment was lower at some sites, but the sample sizes were small.

E. Pivotal Study for JUVÉDERM[®] VOLUMA[®] XC for Correction of Temple Hollowing

Study Design

A multi-center, evaluator-blinded, randomized, parallel-group, controlled pivotal clinical study was conducted to evaluate the safety and effectiveness of JUVÉDERM® VOLUMA® XC for correction of temple hollowing. Subjects were randomized to treatment or no-treatment (control) groups in a 2:1 ratio. Treatment group subjects were given initial study treatment and an optional touch-up treatment 1 month after initial treatment, if deemed necessary to achieve optimal improvement, with follow-up visits at 1, 3, 6, 9, and 13 months after the last treatment. Maintenance treatment was offered to treatment group subjects at 13 months, with follow-up visits 1, 3, and 6 months after treatment. Control group subjects participated in follow-up visits at 1 and 3 months during the "no treatment" control period. Thereafter, control subjects were offered study treatment and optional touch-up with posttreatment follow-up visits at 1, 3, 6, 9, and 13 months after last treatment.

The TI determined the appropriate volume of JUVÉDERM® VOLUMA® XC to be injected in the supraperiosteal plane of the temple area of all subjects. To determine the injection site, the TI identified the point where the temporal fusion line meets the lateral orbital rim and moved superolaterally approximately 1 cm and then inferolaterally approximately 1 cm (Figure 5). The TI inspected and palpated the injection site to identify the location of any superficial vessels, including the superficial temporal artery. The precise injection site was to be slightly adjusted according to the subject's unique anatomy to lower the risk of injecting directly into vessels. The TI inserted the needle at the injection site at an 85- to 90-degree angle to the plane of the temple until it reached the periosteum. The TI injected a single bolus of JUVÉDERM® VOLUMA® XC into the treatment area slowly using gentle, even pressure on the syringe. After injection, if more than one bolus was needed to achieve optimal results, the needle was repositioned in the injection area by changing the injection angle or moving from the initial injection point. No more than 3 different needle entry points were used for a single temple in a single treatment session.

Figure 5. Temple Area Treated



Study Endpoints

With regard to safety, electronic diaries were used by subjects after treatment to record specific signs and symptoms experienced during each of the first 30 days after initial treatment (daily for the first 14 days and every other day for the final 16 days) and for 14 days after touch-up and repeat treatments. Subjects were instructed to rate each treatment site response listed on the diary as "Mild (easily tolerated)," "Moderate (affecting daily activity)," "Severe (unable to do daily activity)," or "None." AEs were reported by the TI at all follow-up visits.

With regard to effectiveness, the 3-month primary effectiveness measure for the study was based on the blinded El's live assessment of temple hollowing using the validated 5-point photonumeric Allergan Temple Hollowing Scale (ATHS) (Table 32). A responder was defined as a subject with \geq 1-point improvement in the ATHS score from baseline. Effectiveness of JUVÉDERM® VOLUMA® XC was demonstrated if at least 60% of subjects treated with JUVÉDERM® VOLUMA® XC were responders at Month 3, and if the responder rate for the treatment group was statistically superior to that of the notreatment control group at Month 3.

Table 32: Allergan Temple Hollowing Scale

Score	Grade	Description
0	Convex	Rounded temple
1	Flat	Flat temple; temporal fusion line may be visible
2	Minimal	Shallow depression or concavity with minimal volume loss; temporal fusion line may be visible
3	Moderate	Moderate depression or concavity with moderate volume loss; moderate prominence of temporal fusion line
4	Severe	Deeply recessed, sunken appearance; marked prominence of temporal fusion line and zygomatic arch

Secondary effectiveness measures included independent, noncollaborative assessments by both the El and the subject of the improvement in the temples using the 5-point GAIS. Secondary measures also included the overall score on the *Satisfaction with Facial Appearance* module and the overall score on the *Satisfaction with Facial Appearance* module and the overall score on the *Satisfaction with Temples* module of the validated FACE-Q questionnaire (0 to 100, where higher scores reflect a better outcome) as assessed by the subjects. Effectiveness was demonstrated if the mean score of each module was statistically superior at Month 3 compared to baseline in the treatment group.

Other effectiveness assessments included measurements of temple volume change from 3D images, subject self-perception of age, and subject assessments of satisfaction with the treatment result and natural look and feel, treatment meeting expectations, likelihood of continuing treatment, and assessment of willingness to recommend treatment.

Subject Demographics

Of 205 subjects enrolled in the study, 34 were screen failures primarily due to ineligibility, and 171 were randomized per protocol, with 113 in the treatment group and 58 in the control group (Table 33). Of the 171 randomized subjects, 162 (94.7%; 106 treatment and 56 control) completed the Month 3 primary endpoint visit, and 53 of the 58 control group subjects (91.4%) opted to receive study treatment after the completion of the 3-month control period. A total of 137 (80.1%; 91 treatment and 46 control) completed the follow-up period through 13 months after treatment.

Table 33: Subject Disposition

Dispecition	Numbe	Number of Subjects		
Disposition	Treatment	Control	Total	
Enrolled	N/A	N/A	205	
Screen Failures	N/A	N/A	34	
Randomized Subjects	113	58	171	
Completed Control Period (Month 3 Primary Endpoint)	106	56	162	
Continued After Control Period	106	53	159	
Lost to Follow-up	9	4	13	
Withdrawal by Subject	6	3	9	
Completed Follow-up Period Through 13 Months After Treatment	91	46	137	
Treatment Group – Did Not Receive Optional Maintenance Treatment (Completed Study)	51	N/A	51	
Treatment Group – Received Optional Maintenance Treatment (VMT Population)	40	N/A	40	
Lost to Follow-up	3	N/A	3	
Withdrawal by Subject	2	N/A	2	
Treatment Group – Completed Follow-up Period Through 6 Months After Maintenance Treatment (Completed Study)	35	N/A	35	

The demographics of the study population are typical for a study performed in the US. Subject demographics and pre-treatment characteristics are presented in Table 34.

 Table 34: Subject Demographics and Pre-Treatment Characteristics

	JUVÉDERM®						
	VOLUMA® XC	Control					
	(N = 113) % (n/N)	(N = 58) % (n/N)					
Gender							
Female	85.0% (96/113)	81.0% (47/58)					
Male	15.0% (17/113)	19.0% (11/58)					
Age (years)		•					
Median	54.0	55.5					
Range (min, max)	32, 82	25, 75					
Race							
White	83.2% (94/113)	82.8% (48/58)					
Black or African American	10.6% (12/113)	10.3% (6/58)					
Asian	1.8% (2/113)	1.7% (1/58)					
American Indian or Alaska Native	0.9% (1/113)	0% (0/58)					
Multiple	3.5% (4/113)	5.2% (3/58)					
Ethnicity							
Hispanic or Latino	23.9% (27/113)	24.1% (14/58)					
Not Hispanic or Latino	76.1% (86/113)	75.9% (44/58)					
Fitzpatrick Skin Type							
1/11	34.5% (39/113)	34.5% (20/58)					
III/IV	53.1% (60/113)	58.6% (34/58)					
V/VI	12.4% (14/113)	6.9% (4/58)					
Baseline ATHS Score							
Minimal	16.8% (19/113)	12.1% (7/58)					
Moderate	47.8% (54/113)	56.9% (33/58)					
Severe	35.4% (40/113)	31.0% (18/58)					

Treatment Characteristics

Of the 113 treatment group subjects, 112 subjects received initial treatment, 73 subjects received touch-up treatment, and 40 subjects received maintenance treatment. Injections were administered using 27 G ½" or 25 G 1" needles for supraperiosteal injection of JUVÉDERM® VOLUMA® XC to the temples. The median total injection volume for initial and touch-up treatments combined in the treatment group was 3.65 mL (range, 0.4 to 8.0 mL) for both temples, with a median of 2.0 mL for initial treatment and 1.5 mL for touch-up

treatment. In the control group, 53 subjects received initial treatment and 38 received touch-up treatment, with similar injection volumes to those in the treatment group. In the maintenance treatment group, the median total volume injected for initial and touch-up treatments combined was 3.81 mL (range, 1.1 to 6.0 mL) for both temples and 2.0 mL (range, 0.2 to 4.0 mL) for maintenance treatment.

Primary Effectiveness Results

JUVÉDERM® VOLUMA® XC provided a clinically and statistically significant correction of moderate to severe temple hollowing compared to the no-treatment control group. The primary effectiveness endpoint was met as (1) greater than 60% of subjects in the treatment group were responders (80.4% improved by \geq 1 point compared with their pre-treatment assessment), and (2) the responder rate for the treatment group was significantly greater (*p* < 0.0001) than the responder rate for the control group (a difference of 66.9%) at Month 3 (Table 35). JUVÉDERM® VOLUMA® XC was found to be effective in all Fitzpatrick skin phototype subgroups, for males and females, and across the studied age range.

Table 35: Effectiveness Summary Responder Rate at 3 Months Based on Evaluating Investigators' Assessments

	Responder Rate at Month 3	p-value
Treatment Group	80.4% (90.8*/113)	N/A
Control Group	13.5% (7.8*/58)	N/A
Difference in Responder Rates (Treatment rate - Control rate)	66.9%	< 0.0001

* Average across datasets with multiple imputation for missing data.

JUVÉDERM[®] VOLUMA[®] XC provided immediate improvement in moderate to severe temple hollowing that lasted for more than 1 year. The ATHS responder rate for the treatment group was 81.7% (85/104) at the 1-month follow-up visit and 73.3% (66/90) at the 13-month follow-up visit.

Secondary Effectiveness Results

The blinded El and subject GAIS responder rates at Month 3 were 83.8% (88/105) and 92.9% (92/99), respectively, for the treatment group. The El and subject GAIS responder rates were 10.7% (6/56) and 2.0% (1/51), respectively, for the untreated control group. The El and subject GAIS responder rates remained high for the treatment group from Month 1 (86.5% and 91.4%, respectively) through Month 13 (72.2% and 78.9%, respectively).

Subject satisfaction with facial appearance, as measured by the FACE-Q *Satisfaction with Facial Appearance* questionnaire, demonstrated significant improvements after treatment with overall mean score of 41.5 at baseline improving to a mean score of 66.3 at Month 3 (p < 0.0001). In the no-treatment control group, the mean score was 39.6 at baseline and decreased to a mean score of 32.7 at Month 3. Improvements were seen from baseline to Month 3 in treatment group responses of *somewhat* or *very satisfied* in the individual FACE-Q *Satisfaction with Facial Appearance* questions shown below:

- How symmetric their face looks: 54.0% of subjects were somewhat or very satisfied at baseline, which improved to 89.0% at Month 3
- How balanced their face looks: 49.6% improved to 88.0%
- How well-proportioned their face looks: 54.0% improved to 87.0%
- How their face looks at the end of the day: 33.6% improved to 72.0%
- How fresh their face looks: 23.0% improved to 68.0%
- How rested their face looks: 25.0% improved to 72.7%
- How their face looks in profile (side view): 31.3% improved to 71.7%
- How their face looks in photos: 31.3% improved to 74.7%
- How their face looks when they first wake up: 29.5% improved to 67.7%
- How their face looks under bright lights: 23.2% improved to 64.6%

At Month 13, the majority of treatment group subjects (more than 62%) were still *somewhat* or *very satisfied* on each of the above questions.

Subject satisfaction with temples, as measured by the FACE-Q *Satisfaction with Temples* questionnaire, also demonstrated significant improvements after treatment with overall mean score of 34.5 at baseline improving to a mean score of 77.4 at Month 3 (p < 0.0001). In the no-treatment control group, the mean score was 30.9 at baseline and decreased to a mean score of 26.5 at Month 3. Improvements were seen from baseline to Month 3 in treatment group responses of *somewhat* or *very satisfied* in the individual FACE-Q *Satisfaction with Temples* questions shown below:

- How their temples look compared with other people their age: 33.9% of subjects were somewhat or very satisfied at baseline, which improved to 86.7% at Month 3
- How well their temples fit in with the rest of their face (balanced): 28.6% improved to 85.7%
- How well the shape of their temples compliments the shape of their face: 23.2% improved to 84.7%
- How their temples look in photos: 33.0% improved to 89.9%
- How their temples look under a bright light: 24.1% improved to 86.9%
- The age their temples make them look: 17.9% improved to 84.7%
- How their temples look when they turn their face slightly to the side: 21.4% improved to 89.9%
- How their temples look from the side (profile view): 29.2% improved to 88.9%
- The shape of their temples: 22.1% improved to 87.0%
- How their temples look in a mirror (straight on): 23.0% improved to 89.0%
- How full their temples look: 17.7% improved to 86.0%
- How youthful their temples make them look: 13.4% improved to 83.7%

At Month 13, the majority of treatment group subjects were still *somewhat* or *very satisfied* on each of the above questions.

Other Effectiveness Results

Other effectiveness assessments included measurements of temple volume change from 3D images, subject self-perception of age, and subject assessments of satisfaction with the treatment result, natural look and feel, treatment meeting expectations, likelihood of continuing treatment, and assessment of willingness to recommend treatment. These effectiveness assessments showed high levels of satisfaction based on subject-reported improvements, and objective measurements showed an increase in volume compared to baseline (Table 36). Improvement continued through 13 months.

Table 36: Additional Treatment Group Effectiveness Results

Assessment	Month 3	Conclusion
Temple Volume Change (Mean Change per Temple)	1.0 cc	The measurement showed an increase after treatment
Subject Self-Perception of Age (Mean Years Younger Than Baseline)	5 years	Subjects reported looking younger after treatment
Subject Satisfaction with Treatment (% Satisfied or Definitely Satisfied)	91.0%	Subjects reported satisfaction with treatment
Subject Satisfaction with Natural Look (% Satisfied or Definitely Satisfied)	89.0%	Subjects reported satisfaction with natural look after treatment
Subject Satisfaction with Natural Feel (% Satisfied or Definitely Satisfied)	91.0%	Subjects reported satisfaction with natural feel after treatment
Subject Assessment of Treatment Meeting Expectations (% responding yes)	87.0%	Subjects reported treatment met expectations
Subject Likelihood of Continuing Treatment (% responding yes)	93.0%	Subjects reported likely to continue treatment
Subject Willingness to Recommend Treatment (% responding yes)	98.0%	Subjects reported willingness to recommend treatment

Follow-up After Maintenance Treatment

Maintenance treatment with JUVÉDERM® VOLUMA® XC was administered to 40 subjects. The effectiveness profile after maintenance treatment was similar to that after initial treatment. The ATHS responder rate after maintenance treatment was 82.9% (29/35) at 1 month, 73.5% (25/34) at 3 months, and 71.4% (25/35) at 6 months.

Subgroup Analysis

The baseline ATHS severity subgroup was evaluated for potential association with outcomes and the effectiveness and safety data are presented in Table 37 and Table 38, respectively. JUVÉDERM® VOLUMA® XC was found to be safe and effective in all Fitzpatrick skin phototype subgroups, for males and females. Additionally, JUVÉDERM® VOLUMA® XC was evaluated for safety and effectiveness across the subject age subgroups (Table 39, Table 40, and Table 41). Although the responder rate for the 30-39 age group is lower than for the rest of the populations, it is difficult to draw any conclusions because of the small population in this subgroup.

		Baseline ATHS Subgroup				
Assessment	Group	Minimal	Moderate	Severe		
ATHS Responder Rate,	JUVÉDERM® Voluma® XC	9/17 (52.9%)	43/51 (84.3%)	34/37 (91.9%)		
n/N (%)	Control	1/6 (16.7%)	3/32 (9.4%)	3/18 (16.7%)		
Investigator GAIS Responder Rate, n/N (%)	JUVÉDERM® Voluma® XC	12/17 (70.6%)	45/51 (88.2%)	31/37 (83.8%)		
nesponder nate, II/N (76)	Control	2/6 (33.3%)	2/32 (6.3%)	2/18 (11.1%)		
Subject GAIS Responder	JUVÉDERM® Voluma® XC	15/16 (93.8%)	44/48 (91.7%)	33/35 (94.3%)		
Rate, n/N (%)	Control	0/4 (0.0%)	0/30 (0.0%)	1/17 (5.9%)		
FACE-Q Satisfaction With Facial Appearance Mean	JUVÉDERM® Voluma® XC	26.4 (17)	21.5 (48)	26.7 (35)		
Change From Baseline (n)	Control	-7.5 (4)	-7.2 (30)	-5.2 (17)		
FACE-Q Satisfaction With Temple Mean Change	JUVÉDERM® Voluma® XC	42.9 (17)	40.9 (48)	45.0 (35)		
From Baseline (n)	Control	-4.0 (4)	-5.1 (30)	-2.6 (17)		

Table 38: Safety Results by Baseline ATHS Subgroups (JUVÉDERM® VOLUMA® XC Treated Population)

Assessment	Baseline ATHS Subgroup						
ASSESSMENT	Minimal	Moderate	Severe				
Any Treatment-related TEAEs, n/N (%)	5/24 (20.8%)	16/84 (19.0%)	8/57 (14.0%)				
Injection Site Discomfort	0/24 (0.0%)	2/84 (2.4%)	0/57 (0.0%)				
Injection Site Pain	1/24 (4.2%)	1/84 (1.2%)	1/57 (1.8%)				
Injection Site Joint Movement Impairment	0/24 (0.0%)	1/84 (1.2%)	0/57 (0.0%)				
Injection Site Mass	0/24 (0.0%)	1/84 (1.2%)	1/57 (1.8%)				
Injection Site Paraesthesia	0/24 (0.0%)	1/84 (1.2%)	0/57 (0.0%)				
Implant Site Pain	0/24 (0.0%)	0/84 (0.0%)	1/57 (1.8%)				
Injection Site Joint Discomfort	0/24 (0.0%)	0/84 (0.0%)	1/57 (1.8%)				
Injection Site Oedema	0/24 (0.0%)	0/84 (0.0%)	1/57 (1.8%)				
Pain in Jaw	4/24 (16.7%)	4/84 (4.8%)	2/57 (3.5%)				
Temporomandibular Joint Syndrome	0/24 (0.0%)	2/84 (2.4%)	0/57 (0.0%)				
Trismus	0/24 (0.0%)	2/84 (2.4%)	0/57 (0.0%)				
Mastication Disorder	0/24 (0.0%)	1/84 (1.2%)	0/57 (0.0%)				
Headache	0/24 (0.0%)	4/84 (4.8%)	4/57 (7.0%)				
Vasodilatation	0/24 (0.0%)	0/84 (0.0%)	1/57 (1.8%)				
Any Serious Adverse Events (SAE), n/N (%)	0/24 (0.0%)	4/84 (4.8%)	3/57 (5.3%)				
Any Treatment Site Responses (TSRs) after Initial Treatment, n/N (%)	19/24 (79.2%)	42/80 (52.5%)	34/57 (59.6%)				

Table 39: Effectiveness Results at Month 3 by Subject Age Subgroups (Initial Treatment Population in Control Period)

Effectiveness			Subject Ag	je Subgrou	ps (mITT F	opulation)	
Measures	Group	30 – 39 years			60 – 69 years	70 – 79 years	≥ 80 years
ATHS Responder	JUVÉDERM® Voluma® XC	50.0% (4/8)	81.0% (17/21)	85.4% (35/41)	88.5% (23/26)	75.0% (6/8)	100.0% (1/1)
Rate, % (n/N)	Control	0.0% (0/4)	30.0% (3/10)	3.8% (1/26)	15.4% (2/13)	33.3% (1/3)	NA
Investigator GAIS	JUVÉDERM® Voluma® XC	75.0% (6/8)	90.5% (19/21)	78.0% (32/41)	84.6% (22/26)	100.0% (8/8)	100.0% (1/1)
Responder Rate, % (n/N)	Control	25.0% (1/4)	30.0% (3/10)	3.8% (1/26)	7.7% (1/13)	0.0% (0/3)	NA
Subject GAIS Responder Rate,	JUVÉDERM® Voluma® XC	100.0% (8/8)	95.2% (20/21)	90.0% (36/40)	91.7% (22/24)	100.0% (5/5)	100.0% (1/1)
% (n/N)	Control	0.0% (0/4)	10.0% (1/10)	0.0% (0/22)	0.0% (0/12)	0.0% (0/3)	NA
FACE-Q Satisfaction With Facial	JUVÉDERM® Voluma® XC	16.5 (8)	16.2 (21)	29.7 (41)	26.3 (24)	8.0 (5)	56.0 (1)
Appearance Mean Change From Baseline (n)	Control	-17.3 (4)	-5.8 (10)	-6.4 (22)	-3.8 (12)	-7.3 (3)	NA
FACE-Q Satisfaction With	JUVÉDERM® Voluma® XC	31.0 (8)	38.0 (21)	47.9 (41)	41.9 (24)	40.2 (5)	51.0 (1)
Temple Mean Change From Baseline (n)	Control	-11.5 (4)	-8.4 (10)	-1.9 (22)	0.9 (12)	-17.3 (3)	NA

Table 40: Treatment Site Responses After Initial Treatment by Subject Age Subgroups (JUVÉDERM® VOLUMA® XC Treated Population)

Age Subgroups	Any TSR %	Pain After Injection %	Tenderness to Touch %	Redness %			Bruising %	Discoloration %	Itching %	
Total (n = 161)	59.0%	50.9%	49.7%	41.0%	37.9%	30.4%	26.1%	16.1%	6.8%	3.7%
30 – 39 years (n = 11)	81.8%	63.6%	72.7%	45.5%	45.5%	18.2%	18.2%	9.1%	0.0%	0.0%
40 – 49 years (n = 30)	56.7%	53.3%	50.0%	26.7%	43.3%	36.7%	33.3%	16.7%	3.3%	3.3%
50 – 59 years (n = 69)	58.0%	53.6%	49.3%	47.8%	31.9%	27.5%	27.5%	17.4%	10.1%	5.8%
60 - 69 years (n = 40)	62.5%	50.0%	50.0%	42.5%	42.5%	40.0%	25.0%	20.0%	7.5%	2.5%
70 – 79 years (n = 10)	30.0%	10.0%	20.0%	20.0%	30.0%	10.0%	0.0%	0.0%	0.0%	0.0%
\geq 80 years (n = 1)	100.0%	100.0%	100.0%	100.0%	100.0%	0.0%	100.0%	0.0%	0.0%	0.0%

Table 41: Treatment-Related AEs After Initial/Touch-up Treatment by Subject Age Subgroups (JUVÉDERM® VOLUMA® XC Treated Population)

Age Subgroups	Any Treatment-Related AE %	Injection Site Discomfort %	Injection Site Pain %	Injection Site Joint Movement Impairment %	Injection Site Mass %	Injection Site Paraesthesia %	Implant Site Pain%	Injection Site Joint Discomfort %	Injection Site Oedema %	Pain in Jaw %	Temporomandibular Joint Syndrome %	Trismus %	Mastication Disorder %	Headache %	Vasodilatation %
30 – 39 years (n = 12)	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
40 – 49 years (n = 30)	13.3%	3.3%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	10.0%	0.0%	0.0%	0.0%	3.3%	0.0%
50 – 59 years (n = 71)	25.4%	1.4%	4.2%	1.4%	2.8%	1.4%	0.0%	0.0%	0.0%	7.0%	1.4%	1.4%	1.4%	0.0%	0.0%
60 - 69 years (n = 40)	17.5%	0.0%	0.0%	0.0%	0.0%	0.0%	2.5%	2.5%	2.5%	5.0%	2.5%	2.5%	0.0%	5.0%	2.5%
70 – 79 years (n = 11)	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
\geq 80 years (n = 1)	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%

8. INSTRUCTIONS FOR USE

A. To Attach Needle to Syringe

STEP 1: Remove tip cap Hold syringe and pull tip cap off the syringe, as shown in Figure A.



STEP 2: Insert needle

Hold the syringe body and firmly insert the hub of the needle (provided in the JUVÉDERM® VOLUMA® XC package) into the LUER-LOK® end of the syringe.

STEP 3: Tighten the needle

Tighten the needle by turning it firmly in a clockwise direction (see Figure B) until it is seated in the proper position, as shown in Figure C.

NOTE: If the position of the needle cap is as shown in Figure D, it is not attached correctly. Continue to tighten until the needle is seated in the proper position.



Figure B

Figure D X

STEP 4: Remove the needle cap

Hold the syringe body in one hand and the needle cap in the other. Without twisting, pull in opposite directions to remove the needle cap, as shown in Figure E.



B. Health Care Professional Instructions

- JUVÉDERM[®] VOLUMA[®] XC injectable gel is a crosslinked, robust, injectable gel formulation, injected using a 27G ½[°] or 25G 1[°] needle or a 25G 1 ½[°] cannula to volumize and contour the cheek for correction of mid-face volume deficit and to augment the chin region to improve the chin profile.
- JUVÉDERM[®] VOLUMA[®] XC injectable gel is injected using a 27G ¹/₂" or 25G 1" needle to volumize the temple region for correction of temple hollowing.
- 3. The TSK STERIGLIDE® 25G 1 ½" cannula was used in the clinical trials (cannula study for cheek and chin study) and is recommended for use with JUVÉDERM® VOLUMA® XC. An entry point was made in the skin with the TSK 23 G introducer needle. In the chin clinical study, JUVÉDERM® VOLUMA® XC was injected into the pogonion (needle), menton (needle or cannula), and left and right pre-jowl sulci (needle or cannula). However, the safety of JUVÉDERM® VOLUMA® XC with cannula for cheek augmentation has not been established in Fitzpatrick Skin Types V and VI. The safety of JUVÉDERM® VOLUMA® XC with cannula for the correction of temple hollowing has not been studied.
- 4. JUVÉDERM® VOLUMA® XC with needle was studied in all Fitzpatrick Skin Types for deep (subcutaneous and/or supraperiosteal) injection for cheek augmentation to correct age-related volume deficit in the mid-face and supraperiosteal injection for correction of moderate to severe temple hollowing in adults over the age of 21.
- 5. Educational resources are available through Allergan Aesthetics, including training in facial anatomy and vasculature, safe injection techniques, and identification and management of potential adverse events, including vascular complications. Health care professionals may contact Allergan Aesthetics for educational and training resources.
- Completion of device-specific use training is required and will be verified prior to distribution of JUVÉDERM[®] VOLUMA[®] XC to improve temple hollowing.
- 7. Before and after treatment, health care professionals are encouraged to conduct vision assessments, including visual acuity, extraocular motility, and visual field testing.
- 8. Health care practitioners are encouraged to be prepared with the following in the event of an intravascular injection:
 - Ensuring supplies are immediately available, as recommended by the American Society for Dermatologic Surgery guidelines¹
 - Identifying a local ophthalmologist or ophthalmology subspecialist to be available in the event of an ophthalmic adverse event related to a dermal filler injection
 - Conducting a basic neurologic examination in the event of an ophthalmic adverse event due to the association of such events with central nervous system deficits
- 9. Prior to treatment, the patient's medical history should be obtained, and the patient should be fully apprised of the indications, contraindications, warnings, precautions, treatment responses, adverse reactions, and method of administration. Patients also should be advised that supplemental "touch-up" implantations may be required to achieve and maintain maximum correction.
- 10. The patient's soft-tissue deficiencies should be characterized with regard to etiology, distensibility, stress at the site, and depth of lesion. Pre-treatment photographs are recommended.
- 11. Topical or injectable anesthesia may be used to manage pain during and after injection.
- 12. After ensuring that the patient has thoroughly washed the treatment area with soap and water, the area should be prepped with alcohol or other antiseptic. Prior to injecting, depress the plunger rod until the product flows out of the needle/cannula.
- 13. If the needle/cannula is blocked, do not increase the pressure on the plunger rod. Instead, stop the injection and replace the needle/cannula.

- 14. When using a cannula, an entry point is made in the skin, e.g. with a sharp needle of appropriate size.
- 15. After insertion of the needle, and just before injection, the plunger rod should be withdrawn slightly to aspirate and verify the needle is not intravascular.
- 16. After the first small amount of material has been injected into the patient, wait a full 3 seconds to allow the lidocaine to take effect before proceeding with the rest of the injection.
- 17. The injection technique for JUVÉDERM® VOLUMA® XC with regard to the angle and orientation of the bevel, the depth (subcutaneous and/or submuscular/supraperiosteal) of injection, and the quantity administered may vary depending on the area being treated. Injection of JUVÉDERM® VOLUMA® XC too superficially (intradermally), or in large volumes over a small area, may result in visible and persistent lumps and/or discoloration.
- 18. Tunneling, fanning, crosshatching, and ferning techniques may be used with a needle or cannula to deliver JUVÉDERM® VOLUMA® XC to achieve optimal results. Serial puncture may be used with a needle to deliver JUVÉDERM® VOLUMA® XC to achieve optimal results. Bolus and serial puncture techniques may be used to achieve optimal results in the chin area. Injection may be administered in an antegrade or retrograde fashion. Inject JUVÉDERM® VOLUMA® XC while applying even pressure on the plunger rod and slowly moving the needle/cannula in the subcutaneous or submuscular/supraperiosteal plane.
- 19. For correction of mid-face volume deficit and augmentation of the chin region, JUVÉDERM® VOLUMA® XC should be distributed in small aliquots (small boluses of 0.1 mL to 0.2 mL) over a large area to reduce the risk of persistent lumpiness.
- 20. For improvement of temple hollowing, JUVÉDERM[®] VOLUMA[®] XC should be administered as a bolus (up to a maximum of 3 boluses per side during a single visit) to the supraperiosteal plane.
- 21. With submuscular/supraperiosteal injection, the number of times the needle passes through the muscle should be minimized to reduce the risk of bruising. It is important to stop injecting before the needle tip reaches the level of the deep dermis to prevent material from being placed too superficially in the skin.
- 22. Correct to 100% of the desired volume effect. Do not overcorrect. The degree and duration of the correction depend on the character of the defect treated, the tissue stress at the implant site, the depth of the implant in the tissue, and the injection technique. Markedly indurated defects may be difficult to correct.
- 23. If immediate blanching occurs, the injection should be stopped, and the area massaged until it returns to a normal color. Blanching may represent a vessel occlusion. If normal skin coloring does not return, do not continue with the injection. Treat in accordance with American Society for Dermatologic Surgery guidelines, which include hyaluronidase injection.¹
- 24. The area of lost facial volume should be lifted by the end of the injection. When injection is completed, the treated site may be gently massaged to mold the product to the contour of the surrounding tissue and assure that it is evenly distributed and conforms to the contour of the surrounding tissues. If overcorrection occurs, massage the area between your fingers or against an underlying superficial bone to obtain optimal results.
- 25. With patients who have localized swelling, the degree of correction is sometimes difficult to judge at the time of treatment. In these cases, it is better to invite the patient back to the office for a touch-up treatment.
- 26. After the initial treatment, an additional treatment may be necessary to achieve the desired level of correction. The same procedure should be repeated until a satisfactory result is obtained. The need for an additional treatment may vary from patient to patient and is dependent upon a variety of factors such as mid-face volume deficit, chin retrusion, or temple hollowing severity, skin elasticity, and dermal thickness at the treatment site.

- 27. Patients may experience treatment site responses, which typically resolve within 3 days for treatment in the temple and 2 to 4 weeks for treatment in the cheek and the chin. Ice may be applied for a brief period following treatment to minimize swelling and reduce pain.
- The health care professional should instruct the patient to promptly report any evidence of problems possibly associated with the use of JUVÉDERM[®] VOLUMA[®] XC.

C. Patient Instructions

It is recommended that the following information be shared with patients:

- Within the first 24 hours, patients should avoid strenuous exercise and extensive sun or heat exposure. Exposure to any of the above may cause temporary redness, swelling, and/or itching at the treatment sites
- If the treated area is swollen, an ice pack may be applied to the site for a short period
- To report an adverse reaction, phone the Allergan Product Surveillance Department at (877) 345-5372

9. HOW SUPPLIED

JUVÉDERM® VOLUMA® XC injectable gel is supplied in individual treatment syringes with needles as indicated on the carton. JUVÉDERM® VOLUMA® XC can be injected with either a 27G ½¹⁰ needle, a 25G 1¹¹ needle or a 25G 1 ½¹¹ cannula. The TSK *STERiGLIDE*® 25G 1½¹¹ cannula is not supplied with JUVEDÉRM® VOLUMA® XC but is available for purchase through AbbVie. The volume in each syringe is as stated on the syringe label and on the carton. The contents of the syringe are sterile and non-pyrogenic. Do not resterilize. Do not use if package is open or damaged.

10. SHELF LIFE AND STORAGE

JUVÉDERM® VOLUMA® XC injectable gel must be used prior to the expiration date printed on the label. Store at room temperature (up to 25° C/77°F). DO NOT FREEZE.

JUVÉDERM[®] VOLUMA[®] XC injectable gel has a clear appearance. In the event that a syringe contains material that is not clear, do not use the syringe; notify Allergan Product Surveillance immediately at (877) 345-5372.

To place an order, contact Allergan at (800) 377-7790.



Irvine, CA 92612 USA

1-800-624-4261

Made in France

© 2024 AbbVie. All rights reserved.

JUVÉDERM, VOLUMA, and their designs are trademarks of Allergan Holdings France SAS, an AbbVie company, or its affiliates. All other trademarks are the property of their respective owners. Patented. See: www.abbvie.com/patents

atented. See: www.abbvie.com/paten

hcp.Juvederm.com 11/2023 024282 DN7084

¹Alam M, Gladstone H, Kramer EM, et al. ASDS guidelines of care: injectable fillers. *Dermatol Surg.* 2008;34(suppl 1):S115-S148.